



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Flamingo

mm (inch) kg (lb)

	 1	 2	 3	 4
 Overall width, seat	350 mm (13¾")	380 mm (15")	430 mm (17")	480 mm (19")
 Length, seat	260 mm (10¼")	315 mm (12½")	360 mm (14¼")	425 mm (16¾")
 Overall height, seat	400 mm (15¾")	430 mm (17")	480 mm (19")	510 mm (20")
 Effective seat width	240 mm (9½")	280 mm (11")	320 mm (12½")	340 mm (13½")
 Effective seat depth	250 mm (9¾")	290 mm (11½")	350 mm (13¾")	400 mm (15¾")
 Back width, top of backrest	200 mm (7¾")	220 mm (8¾")	260 mm (10¼")	280 mm (11")
 Backrest height	350 mm (13¾")	380 mm (15")	420 mm (16½")	470 mm (18½")
 Folded height, seat	240 mm (9½")	260 mm (10¼")	270 mm (10¾")	270 mm (10¾")
 Folded length, seat	370 mm (15")	400 mm (15¾")	430 mm (17")	460 mm (18")
 Backrest angle	30°	30°	30°	30°
 Mass device, seat	2.5 kg (5.5 lb)	3 kg (6.5 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max rated load	39.5 kg (87 lb)	55 kg (121.2 lb)	66 kg (145.5)	76.5 kg (168.6 lb)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1274.3 Pricing and specifications subject to change without notice.

kg (lb)	1	2	3	4
 Mass device, frame	11,5 kg (25.3 lb)	11,5 kg (25.3 lb)	12,5 kg (27.5 lb)	12,5 kg (27.5 lb)
 Mass device, seat	2,5 kg (5.5 lb)	3 kg (6.6 lb)	4 kg (8.8 lb)	4.5 kg (9.9 lb)
 Mass device	14 kg (30.8 lb)	14,5 kg (31.9 lb)	16,5 kg (36.3 lb)	17 kg (37.4 lb)
 Max user mass	35 kg (77.1 lb)	50 kg (110.2 lb)	60 kg (132.2 lb)	70 kg (154.3 lb)
 Max. load, operating force	35 kg (77.1 lb)	35 kg (77.1 lb)	55 kg (121.2 lb)	55 kg (121.2 lb)
 Max rated load	51 kg (112.4 lb)	66,5 kg (146.6 lb)	78,5 kg (173 lb)	89 kg (196 lb)

Base Frame and Mobility Base Options

Seat and Back



Seat and Back

Seat and Back, Blue, Size 1 (880101-61)	\$754
Seat and Back, Blue, Size 2 (880102-61)	\$806
Seat and Back, Blue, Size 3 (880103-61)	\$858
Seat and Back, Blue, Size 4 (880104-61)	\$910

Flamingo Mobility Base

Footrest only included with the High-low frame



Flamingo Frame W/ 3" Casters, Size 1 (880201-NL)	\$884
Flamingo Frame W/ 3" Casters, Size 2 (880202-NL)	\$931
Flamingo Frame W/ 3" Casters, Size 3 (880203-NL)	\$1014
Flamingo Frame W/ 3" Casters, Size 4 (880204-NL)	\$1087



Potty Frame, Size 1 (880160-1)	\$468
Potty Frame, Size 2 (880160-2)	\$515



High-Low Bath Frame with 4" Casters, Size 1 (880411)	\$2460
High-Low Bath Frame with 4" Casters, Size 2 (880412)	\$2460
High-Low Bath Frame with 4" Casters, Size 3 (880413)	\$2460
High-Low Bath Frame with 4" Casters, Size 4 (880414)	\$2527



Fitting for Direct Toilet Mounting

Direct Toilet Mount, Size 1 (880141)	\$650
Direct Toilet Mount, Size 2 (880142)	\$650
Direct Toilet Mount, Size 3 (880143)	\$650
Direct Toilet Mount, Size 4 (880144)	\$650



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Footrest for Mobility Base

If wanting a footrest for the flamingo frame or Direct Toilet Mounting, please quote the footrests below



Flamingo Foot Support Pair, Size 1 (880211)	\$312
Flamingo Foot Support Pair, Size 2 (880212)	\$338
Flamingo Foot Support Pair, Size 3 (880213)	\$364
Flamingo Foot Support Pair, Size 4 (880214)	\$390



Calf and Footrest Supports

Calf Support (880199) _____ \$74

Foot Straps, Plastic Lock (Pair) (85096) (K0040) _____ \$42



Pelvic and Trunk Positioning

Pelvic Positioning



Seat with Upholstery

Upholstery, Seat w/wide sides, Blue, Size 2 (880166-61)	\$118
Upholstery, Seat w/wide sides, Blue, Size 3 (880167-61)	\$123
Upholstery, Seat w/wide sides, Blue, Size 4 (880168-61)	\$128



Seat Infill

Seat Infill for Seat, Blue, Size 1 (880185-61)	\$61
Seat Infill for Seat, Blue, Size 2 (880186-61)	\$61
Seat Infill for Seat, Blue, Size 3 (880187-61)	\$68
Seat Infill for Seat, Blue, Size 4 (880188-61)	\$68

Splashguard

Splashguard, Low Profile, Size 1, Blue (880195-61)	\$85
Splashguard, Low Profile, Size 2,3,4, Blue (880196-61)	\$85
Splashguard, High Profile, Size 1, Blue (880191-61)	\$90
Splashguard, High Profile, Size 2,3,4, Blue (880192-61)	\$90

Accessories

Seat Belt w/ Plastic Lock (85100)	\$34
Suction Cups for Seat (880150)	\$149

Neoprene Hip Belts

Neoprene Hip Belts, size 1 (85423-1)	\$107
Neoprene Hip Belts, size 2 (85423-2)	\$107
H-Harness with Plastic Lock (85105)	\$68

Trunk Positioning



Sides with Upholstery

Sides for Back W/ Upholstery, Blue (Each), Size 1 (880111-61)	\$143
Sides for Back W/ Upholstery, Blue (Each), Size 2 (880112-61)	\$143
Sides for Back W/ Upholstery, Blue (Each), Size 3 (880113-61)	\$143
Sides for Back W/ Upholstery, Blue (Each), Size 4 (880114-61)	\$158



Back Extension

Trunk Positioning, Back Extension, Blue Size 1 and 2 (880118-61)	\$229
Back Extension, Size 3 and 4 (880119-61)	\$229

Back Reinforcement Kit

Back Reinforcement Kit 3 and 4 (880155)	\$284
---	-------

Cross vest, Water-proof mesh (net)

Cross vest, Water-proof mesh, Size 0 (89625-040)	\$218
Cross vest, Water-proof mesh, Size 1 (89625-140)	\$218
Cross vest, Water-proof mesh, Size 2 (89625-240)	\$218.28
Cross vest, Water-proof mesh, Size 3 (89625-340)	\$220
Cross vest, Water-proof mesh, Size 4 (89625-440)	\$220

Mounting Rod

Rod needed for Back Extension

Mounting Rod for Back Extension L:8" (880120)	\$51
---	------

Tilt Support Rod

Tilt Support Rod, Size 3 and 4	\$156
--------------------------------	-------

Vest, Water-Proof Mesh (Each)

More options can be found at etac.com

Vest, Water-Proof Mesh (Each), Size 3/4 (89111-07540)	\$219
Vest, Water-Proof Mesh (Each), Size 1 (89111-140)	\$219
Vest, Water-Proof Mesh (Each), Size 2 (89111-240)	\$219
Vest, Hygiene, Size 3, 2.75" x 4.5" x 11" x 5.75" x 18" (89111-340)	\$219
Vest, Hygiene, Size 4, 2.75" x 6.5" x 13.5" x 7.5" x 19" (89111-440)	\$219



Phone (877) 227-0097 | Fax (704) 882-0751

Flamingo

Fixed Side Supports



Fixed Side Supports, Gray Cover, Set, Size 1, 3.5" x 4" (88401-10)	\$260
Fixed Side Supports, Gray Cover, Set, Size 2, 4.75" x 5.5" (88401-20)	\$260

Swing-Away Brackets

Swing-Away Bracket, Right, Each (83300-R)	\$250
Swing-Away Bracket, Right, Offset (83500-R)	\$250
Swing-Away Bracket, Left, Each (83300-L)	\$250
Swing-Away Bracket, Left, Offset (83500-L)	\$250

Headrests and Hardware

Head Support No. 1, Blue PU (6071000)	\$228
Head Support No. 2, Blue PU (6072000)	\$228
Head Support No. 3, Blue PU (6073000)	\$228
Head Support No. 4, Blue PU (6074000)	\$228

Cushion, Curved

Cushion, Curved, Blue, Each, L: 5.75" x H: 2.5" (6075409)	\$118
Cushion, Curved, Blue, Each, L: 6.5" x H: 2.5" (6075408)	\$118
Cushion, Curved, Blue, Each, L: 7.5" x H: 2.75" (6075407)	\$118

Bar, Straight, Stainless

Bar, Straight, Stainless L: 4" (82569)	N/C
Bar, Straight, Stainless L: 6" (82648)	N/C

Protection Pad Covers

For Swing-Away Brackets

Protection Pad Covers

Protection Pad Covers, Gray (89151)	\$112
-------------------------------------	-------

Armrest Options

Armrests



Swing-Away Armrest/Tray, Blue, Size 1 (880131-61)	\$255
Swing-Away Armrest/Tray, Blue, Size 2 (880132-61)	\$255
Swing-Away Armrest/Tray, Blue, Size 3 (880133-61)	\$255
Swing-Away Armrest/Tray, Blue, Size 4 (880134-61)	\$255



Separate Armrests, Set, 9" (880152)	\$198
-------------------------------------	-------

Accessories

Push Handles



Push Brace (One Size) (880110)	\$107
--------------------------------	-------



For use with High-Low Base

Helper Lift Handle, Size 1 and 2 (880110-HL)	\$95
Helper Lift Handle, Size 3 and 4 (880110-HLL)	\$95

Potty Options



Potty Ring W/ Holders, Size 1 (880121-PR)	\$143
Potty Ring W/ Holders, Size 2 (880122-PR)	\$143
Potty Ring W/ Holders, Size 3 (880123-PR)	\$143
Potty Ring W/ Holders, Size 4 (880124-PR)	\$143



Commode pan w/holder

Commode pan w/holder, Size 1 (880121)	\$133
Commode pan w/holder, Size 2 (880122)	\$133
Commode pan w/holder, Size 3 (880123)	\$133
Commode pan w/holder, Size 4 (880124)	\$133