



Phone (877) 227-0097 | Fax (704) 882-0751

Manatee

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Manatee

mm (inch)	kg (lbs)				
		0	1	2	3
 Effective seat width		330 mm (13")	330 mm (13")	330 mm (13")	400 mm (15¾")
 Effective seat depth		220-250 mm (8¾-9¾")	300-370 mm (11¾-14½")	390-460 mm (15¼-18")	400-470 mm (15¾-18½")
 Backrest height		430 mm (17")	430 mm (17")	580 mm (22¾")	640 mm (25¼")
 Lower leg length		280 mm (11")	280 mm (11")	350 mm (13¾")	350 mm (13¾")
 Seat surface height at front edge		140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")
 Overall width		420 mm (16½")	420 mm (16½")	420 mm (16½")	490 mm (19¼")
 Folded length		1030-1100 mm (40½-43¾")	1100-1170 mm (43-46")	1360-1430 mm (53½-56¾")	1440-1510 mm (56¾-59½")
 Folded height		130 mm (5")	130 mm (5")	130 mm (5")	130 mm (5")
 Mass device		5.5 kg (12.1 lb)	5.5 kg (12.1 lb)	6 kg (13.2 lb)	7 kg (15.4 lb)
 Maximum user mass		50 kg (110.2 lb)	50 kg (110.2 lb)	75 kg (165.3 lb)	75 kg (165.3 lb)
 Maximum rated load		57.5 kg (126.7 lb)	57.5 kg (126.7 lb)	83 kg (182.9 lb)	84 kg (185.1 lb)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1259.2 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751

Manatee

Manatee Seat and Frame Options

Size and Color

Shower seat only



Blue

Size 0, Blue (882005)	\$780
Size 1, Blue (882015)	\$780
Size 2, Blue (882025)	\$780
Size 3, Blue (882035)	\$852

Red

Size 1, Red (882014)	\$780
Size 2, Red (882024)	\$780
Size 3, Red (882034)	\$852

Manatee Frame Options



Rolling Shower Frame

Rolling Shower Frame (882100)	\$2499
-------------------------------	--------



Tub Slider Frame

Tub Slider without Manatee W/ Casters, Size 1 and 2 (MC6000PED)	\$4100
Tub Slider without Manatee W/ Casters, Size 3 (MC6200PED)	\$5100



Basic Manatee Frame

Basic Manatee Frame (882101)	\$390
------------------------------	-------

Accessories and Straps

Accessories



Leg Straps

Leg Straps, Size 0,1 and 2, Blue (883055)	\$86
Leg Straps, Size 0,1 and 2, Red (883054)	\$86
Leg Straps, Size 3, Blue (883105)	\$109
Leg Straps, Size 3, Red (883104)	\$109



Head/ Lateral Support

Head Lateral Support, Size 0,1 and 2, Blue (883015)	\$123
Head Lateral Support, Size 0,1 and 2, Red (883014)	\$123
Head Lateral Support, Size 3, Blue (883065)	\$123
Head Lateral Support, Size 3, Red (883064)	\$123



Pommel

Pommel, Size 0,1,2, Blue (883045)	\$109
Pommel, Size 0,1,2, Red (883044)	\$109
Pommel, Size 3, Blue (883095)	\$109
Pommel, Size 3, Red (883094)	\$109



Torso Belt

Torso Belt, Size 0,1 and 2, Blue (883025)	\$86
Torso Belt, Size 0,1 and 2, Red (883024)	\$86
Torso Belt, Size 3, Blue (883075)	\$86
Torso Belt, Size 3, Red (883074)	\$86