

Company:

Account Number:

Quote:

ATP:

Marked For:

Associated Quote/Order#:

Ship To:

Address:

City:

Zip:

Email/Fax:

Date:

State:

Phone:

Effective: 5/7/2026

Meerkat

Meerkat measurements

	1	2	3
 User height: floor to axil	15¾ – 29½" (400 - 750 mm)	28¾ – 39¾ (720 - 1000 mm)	33 – 51¼" (840 - 1300 mm)
 User height	27½ – 41¼" (700 - 1050 mm)	37½ – 53¼" (950 - 1350 mm)	46 – 70¼" (1170 - 1800 mm)
 Max user mass	44 lb (20 kg)	143.3 lb (65 kg)	198.4 lb (90 kg)
 Length	23½" (595 mm)	31½" (800 mm)	31½" (800 mm)
 Overall width	20¾" (525 mm)	30¾" (780 mm)	30¾" (780 mm)
 Overall height	14¾ – 29½" (375 - 750 mm)	19¾ – 39½" (500 - 1000 mm)	25½ – 51¼" (650 - 1300 mm)
 Mass device	19.8 lb (9 kg)	33.6 lb (15 kg)	34.7lb (15 kg)

Further measurements at our website

Meerkat Frame Option

Frame and Rockers (Forward/Backwards/Sideways)	
Pastel Blue, Size 1 (870011-1RS)	\$1316
Strawberry Pink, Size 1 (870011-2RS)	\$1316
Grey, Size 2 (870022-RS)	\$1448
Grey, Size 3 (870023-RS)	\$1511

Pelvic and Chest Supports

Pelvic and Chest Support

Pelvic Support



Pelvic Support, C: 17.75" - 23.5", Size Small (869240)	\$451
Pelvic Support, C: 23.5" - 29", Size Medium (869241)	\$455
Pelvic Support, C: 27.5" - 37.5", Size Large (869242)	\$472
Pelvic Support, C: 37.5" - 42.5", Size XL (869243)	\$482
Pelvic Support, C: 39.25" - 47.25", Size XXL (869244)	\$490

Chest Support



Chest Support, C: 17.75" - 23.5", Size Small (869480)	\$451
Chest Support, C: 23.5" - 29", Size Medium (869481)	\$455
Chest Support, 27.5" - 37.5", Size Large (869482)	\$472
Chest Support, C: 37.5" - 43.25", Size XL (869483)	\$482
Chest Support, C: 39.25" - 47.25", Size XXL (869484)	\$490



Daily Adjustment Kit, Height

Daily Adjustment Kit, For Chest/Pelvic and Knee (124027)	\$37
--	------

Sandals and Knee Supports

Sandals and Knee Supports



Knee Support

Knee Supports, Size Small (869840)	\$490
Knee Support, Size Medium (869841)	\$505
Knee Supports, Size Large (869842)	\$520



Sandals

Adjustable Sandals for Size 1 Plate, Small (869490)	\$233
Adjustable Sandals for Size 2 Plate, Medium (869491)	\$244
Adjustable Sandals for Size 2 Plate, Large (869492)	\$251



Additional Straps and Fixlocks for Sandals

Straps and Fixlocks, Medium (869497)	\$60
Straps and Fixlocks, Large (869498)	\$63



Daily Adjustment Kit for Knee Support

Daily Adjustment Kit for Knee Support (124026)	\$42
--	------

Accessories

Accessories



Tray, Size 1 (869485)	\$443
Tray, Size 2 (869486)	\$495



Padding for Meerkat Ratchet Straps (117818-PADS) (Pair)	\$56
---	------



Displacement Kit Bracket - Offset Mounting (124025)	\$81
---	------