



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Mustang

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Mustang

mm (inch), kg (lb)

	 1	 2	 3	 4
 Overall width 4" or 8" wheels	520 mm (20½")	580 mm (22¾")	650 mm (25½")	700 mm (27½")
 Inside width	380 mm (15")	450 mm (17¾")	500 mm (19½")	560 mm (22")
 Length with 4" wheels	750 mm (29½")	830 mm (32¾")	-	-
 Length with 8" wheels	820 mm (32¼")	900 mm (35½")	1000 mm (39¼")	1100 mm (43¼")
 Height with 4" wheels	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	-	-
 Height with 8" wheels	510-650 mm (20-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	1030-1250 mm (40½-49¼")
 Seat height with 4" wheels PU seat (accessory)	170-350 mm (6¾-13¾")	240-410 mm (9½-16¼")	-	-
 Seat height with 8" wheels PU seat (accessory)	260-440 mm (10¼-17¼")	330-500 mm (13-19¾")	450-680 mm (17¾-26¾")	620-840 mm (24½-33")
 Seat height with 4" wheels Comfort seat (accessory)	170-430 mm (6¾-17")	240-490 mm (9½-19¼")	-	-
 Seat height with 8" wheels Comfort seat (accessory)	260-520 mm (10¼-20½")	330-580 mm (13-22¾")	450-760 mm (17¾-30")	620-920 mm (24½-36¼")
 Width between handles (accessory)	165-320 mm (6¼-12½")	165-320 mm (6¼-12½")	140-400 mm (6½-15¾")	140-400 mm (6½-15¾")
 Angle, center bar (a1)	30°	30°	25°	25°
 Minimum turning radius	360 mm (14¼")	550 mm (21¾")	600 mm (23½")	650 mm (26½")
 User height, floor to axle	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	870-1300 mm (34¼-51¼")

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1272.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Mustang

Mustang Frame Options

Mustang, Standard Package



Mustang, Standard Package

Mustang, Yellow, 4" Casters, Size 1 (869011)	\$2154
Mustang, Yellow, 8" Casters, Size 1 (869015)	\$2154
Mustang, Yellow, 8" Casters Slim Line Wheels, Size 1 (869015-SL)	\$2154
Mustang, White, 4" Casters, Size 2 (869025-4)	\$2243
Mustang, White, 8" Casters, Size 2 (869024)	\$2243
Mustang, White, 8" Casters Slim Line Wheels, Size 2 (869024-SL)	\$2243
Mustang, Red, 8" Casters, Size 3 (869033)	\$3171
Mustang, Red, 8" Casters Slim Line Wheels, Size 3 (869033-SL)	\$3171
Mustang, Anthracite (Grey), 8" Casters, Size 4 (869042)	\$3628
Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels, Size 4 (869042-SL)	\$3628

Surrounding Supports and Covers

Surrounding Chest Support and Covers

Required Selection



Surrounding Chest Support

Chest Support, 6"- 8.5", Small (869810)	\$520
Chest Support, 9.5" - 12", Medium (869820)	\$595
Chest Support, 11" - 14.75", Large (869830)	\$695

Covers

Cover for Chest Support, Small (869810-160)	\$129
Cover for Chest Support, Medium (869820-260)	\$158
Cover for Chest Support, Large (869830-360)	\$168

Mounting For Chest Support

Low Mounting Kit for Surrounding Chest Support (869800-KIT)	\$260
---	-------

Seat Options

Seats and Slings



Sling Seat

Sling Seat for Small Chest Support, W: 7" x L: 9.5", Size 1 (869570)	\$195
Sling Seat for Medium Chest Support, W: 7" x L: 12.75", Size 2 (869571)	\$195



Mustang Seat

Mustang Seat, Small, Size 1 and 2 (869110)	\$295
Mustang Seat, Large, Size 3 and 4 (869120)	\$325
Mustang Comfort Seat, Small, Size 1 and 2 (869115)	\$475
Mustang Comfort Seat, Large, Size 3 and 4 (869125)	\$495

Additional Support Options

Supports



Pelvic Supports

Mounts on center rod

Surrounding Pelvic Supports, Small (869610)	\$525
Surrounding Pelvic Supports, Medium (869620)	\$565
Surrounding Pelvic Supports, Large (869630)	\$595

Extra Pad

Extra Pad for Chest/Surrounding Pelvic Supports (869225)	\$69
--	------

Sacrum Support

Requires Mustang Seat

Sacrum Support, Small (869250)	\$185
Sacrum Support, Large (869255)	\$275

Hip Supports

Requires Mustang Seat

Hip Support on Straight Bar, Size 1 and 2 (869221)	\$429
Hip Support on Straight Bar, Size 3 and 4 (869231)	\$429
Hip Support on Curved Bar, Size 1 and 2 (869222)	\$429
Hip Support on Curved Bar, Size 3 and 4 (869232)	\$429



Abdominal Support

Abdominal Support (869260)	\$204
----------------------------	-------



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Mustang

Armrest Options

Armrests

Armrests

Armrest for Anterior Use W/ Handles, Small (869350)	\$535
Armrest for Anterior Use W/ Handles, Large (869355)	\$563

Arm Prompts

Arm Prompts for Anterior Use, Small (869310-1)	\$595
Arm Prompts for Anterior Use, Medium (869320-1)	\$595
Arm Prompts for Anterior Use, Large (869325-1)	\$595
Arm Prompts for Posterior Use, Small (869330-1)	\$625
Arm Prompts for Posterior Use, Medium (869340-1)	\$625
Arm Prompts for Posterior Use, Large (869345-1)	\$625



Padded Forearm Positioner

If wanting a pair of straps, quantity two needs to be selected

Padded Forearm Positioner Strap (Each) (869370) _ \$75

Quantity:

Decelerating Handbrake

Decelerating Handbrake, Size 1 (869346-1)	\$1525
Decelerating Handbrake, Size 2 (869346-2)	\$1525
Decelerating Handbrake, Size 3 (869346-3)	\$1525
Decelerating Handbrake, Size 4 (869346-4)	\$1525

Headrest and Hardware

Head Support

Headrest and connecting rod requires the Mustang Seat and either the Sacrum Support or the Hip Supports



Headrest

Anatomic Head Support, Small (5099001)	\$285
Anatomic Head Support, Medium (5099002)	\$285
Anatomic Head Support, Large (5099003)	\$285

Headrest Mounting Hardware

Headrest Mounting Hardware, Size 1 and 2, Short (98619-25)	\$250
Headrest Mounting Hardware, Size 3 and 4, Long (98619-26)	\$256

Headband System

Headband System (99500)	\$273
-------------------------	-------

Headrest Cover

Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium (9945-83)	N/C
Anatomic Head Support Cover, Large (9967-83)	N/C



Connecting Rod

Required for use of headrest

Connecting Rod Between Hip and Head Support (869226) _ \$85

Accessories

Leg Separator and Ankle Prompts



Leg Separators

Leg Separator for Size 1 w/ 4" Wheels (869410)	\$317
Leg Separator for Size 1 w/ 8" Wheels (869410-8)	\$317
Leg Separator for Size 2 w/ 8" Wheels (869420)	\$317
Leg Separator for Size 2 w/ 4" wheels and Size 3 (869430)	\$317
Leg Separator for Size 4 (869440)	\$382



Leg Pad

Leg Pad (869470)	\$147
------------------	-------



Ankle Prompts

Ankle Prompts for Size 1 (869450)	\$225
Ankle Prompts for Size 2 w/ 8" Wheels (869455)	\$225
Ankle Prompts for Size 2 w/ 4" Wheels (869455-4)	\$225
Ankle Prompts for Size 3 (869460)	\$235
Ankle Prompts for Size 4 (869465)	\$235



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Mustang

Accessories



Attendant Guide Bar
Attendant Guide Bar, Size 1 and 2 (869510) _____ \$250



Storage Pouch
Storage Pouch, Size 1 and 2 (869560) _____ \$109



Storage Bag
Storage Bag, 3 and 4 (869565) _____ \$218