



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

x:panda

mm (inch) kg (lb)

	1	2	3	4
Effective seat width	160, 200, 240 mm (6¼", 7¾", 9½")	200, 240, 280 mm (7¾", 9½", 11")	260, 310, 360 mm (10¼", 12¼", 14¼")	340, 400, 460 mm (13½", 15¾", 18")
Width, thighs	70, 90, 110 mm (2¾", 3½", 4¼")	100, 120, 140 mm (4", 4¾", 5½")	120, 145, 170 mm (4¾", 5¾", 6¾")	130, 190, 250 mm (5", 7½", 9¾")
Effective seat depth	160-260 mm (6¼"-10¼")	220-330 mm (8½"-13")	260-390 mm (10¼"-15¼")	390-530 mm (15¼"-20¾")
Backrest height	270-370 mm (10½"-14½")	320-420 mm (12½"-16½")	350-470 mm (13¾"-18½")	500-710 mm (19¾"-28")
Back width, top of backrest	230 mm (9")	280 mm (11")	330 mm (13")	350+50 mm (13¾"+2")
Overall width, both legs in max abduction	310, 350, 390 mm (12¼", 13¾", 15¼")	380, 420, 460 mm (15", 16½", 18")	470, 520, 570 mm (18½", 20½", 22½")	590, 650, 710 mm (23¾", 25½", 28")
Armrest to seat distance	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")
Front location of armrest structure	140 mm (5½")	140 mm (5½")	180 mm (7")	350 mm (13")
Backrest angle	(-5°)-(30°)	(-5°)-(30°)	(-5°)-(30°)	(-5°)-(30°)
Abduction angle, each leg (s1) Size 1 in the smallest position	0°-15° 0°-8°	0°-15°	0°-15°	0°-15°
Maximum user mass	30 kg (66.1 lb)	30 kg (66.1 lb)	60 kg (132.2 lb)	80 kg (176.3 lb)
Maximum user mass at transportation	30 kg (66.1 lb)	30 kg (66.1 lb)	60 kg (132.2 lb)	75 kg (165.3 lb)
Mass device*	6.5 kg (14.3 lb)	8.5 kg (18.7 lb)	11 kg (24.3 lb)	13 kg (28.5 lb)
Max rated load	38.5 kg (84.7 lb)	40.5 kg (89.2 lb)	73 kg (160.9 lb)	95 kg (209.4 lb)
Minimum User height	650 mm (25½")	850 mm (33½")	1000 mm (39½")	1390 mm (54¾")
Maximum User height	1050 mm (41¼")	1280 mm (50½")	1430 mm (56¼")	1800 mm (70¾")

*seat only (without accessories)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1296.3 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

x:panda Base Options

Seat and Back must match Base. x:panda seat comes complete with growth blocks to allow for sm-med-lg setting.

x:panda



Size 1 x:panda

Size 1 x:panda, Small (8910111-02) _____ N/C
Size 1 x:panda, Medium (8910112-02) _____ N/C
Size 1 x:panda, Large (8910113-02) _____ N/C

Size 3 x:panda

Size 3 x:panda, Small (8910131-02) _____ N/C
Size 3 x:panda, Medium (8910132-02) _____ N/C
Size 3 x:panda, Large (8910133-02) _____ N/C

Size 2 x:panda

Size 2 x:panda, Small (8910121-02) _____ N/C
Size 2 x:panda, Medium (8910122-02) _____ N/C
Size 2 x:panda, Large (8910123-02) _____ N/C

Size 4 x:panda

Size 4 x:panda, Small (8910141-02) _____ N/C
Size 4 x:panda, Medium (8910142-02) _____ N/C

x:panda Back and Seat Cushions

x:panda Back

Back selection required. Back size must match above base size.

Size	Price	Small	Medium	Large
Size 1	\$724	Size 1 Back x:panda, Small (X8910111-01) (E2291)	Size 1 Back x:panda, Medium (X8910112-01) (E2291)	Size 1 Back x:panda, Large (X8910113-01) (E2291)
Size 2	\$724	Size 2 Back x:panda, Small (X8910121-01) (E2291)	Size 2 Back x:panda, Medium (X8910122-01) (E2291)	Size 2 Back x:panda, Large (X8910123-01) (E2291)
Size 3	\$724	Size 3 Back x:panda, Small (X8910131-01) (E2291)	Size 3 Back x:panda, Medium (X8910132-01) (E2291)	Size 3 Back x:panda, Large (X8910133-01) (E2291)
Size 4	\$742	Size 4 Back x:panda, Small (X8910141-01) (E2291)	Size 4 Back x:panda, Medium (X8910142-01) (E2291)	Size 4 Back x:panda, Large (X8910143-01) (E2291)

x:panda Seat

Seat selection required. Seat size must match above base size.

Size	Price	Small	Medium	Large
Size 1	\$742	Size 1 Seat x:panda, Small (X8910111-02) (E2292)	Size 1 Seat x:panda, Medium (X8910112-02) (E2292)	Size 1 Seat x:panda, Large (X8910113-02) (E2292)
Size 2	\$742	Size 2 Seat x:panda, Small (X8910121-02) (E2292)	Size 2 Seat x:panda, Medium (X8910122-02) (E2292)	Size 2 Seat x:panda, Large (X8910123-02) (E2292)
Size 3	\$742	Size 3 Seat x:panda, Small (X8910131-02) (E2292)	Size 3 Seat x:panda, Medium (X8910132-02) (E2292)	Size 3 Seat x:panda, Large (X8910133-02) (E2292)
Size 4	\$742	Size 4 Seat x:panda, Small (X8910141-02) (E2292)	Size 4 Seat x:panda, Medium (X8910142-02) (E2292)	

Backrest Options

Must choose a backrest option. Must match base size.

Dynamic Back

Dynamic Back, 100 Newton Gas Spring, Size 1 and 2 (X911100) (K0108) _____ \$784
Dynamic Back, 150 Newton Gas Spring, Size 3 (X8911150) (K0108) _____ \$784
Dynamic Back, 350 Newton Gas Spring, Size 4 (X8911350-4) _____ \$784

Non-Dynamic Back

Non-Dynamic Gas Spring, Size 1 and 2 (8911001) (E1014) _____ \$482
Non-Dynamic Gas Spring, Size 3 (8911002) (E1014) _____ \$482



Recline Option

Required Selection

Recline, '-5° to 30° (X66-8911) (E1014) _____ \$482

Seat Cushion Type

Required selection

Split Seat Cushion _____ N/C

Seat Cushion _____ N/C

For use with Medial Thigh System or Abduction Supports

For use without Medial Thigh System or Abduction Supports



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Back Cushion Upholstery Color

Required selection

Blue Upholstery	N/C
Black Upholstery	N/C
Red Upholstery	N/C
Pink Upholstery	N/C

Not available with size 4

Seat Cushion Upholstery Color

Required Selection

Black Cushion Cover	N/C
Blue Cushion Cover	N/C
Red Cushion Cover	N/C
Pink Cushion Cover	N/C

Not available with size 4

Replacement Black Wipeable Back Cover

Additional Back Cover only. Only available with Split Cushion.

Size	Price	Small	Medium	Large
Size 1	\$128	Back Cover, Black Wipeable, Size 1, Small (8912411-82)	Back Cover, Black Wipeable, Size 1, Medium (8912412-82)	Back Cover, Black Wipeable, Size 1, Large (8912413-82)
Size 2	\$128	Back Cover, Black Wipeable, Size 2, Small (8912421-82)	Back Cover, Black Wipeable, Size 2, Medium (8912422-82)	Back Cover, Black Wipeable, Size 2, Large (8912423-82)
Size 3	\$128	Back Cover, Black Wipeable, Size 3, Small (8912431-82)	Back Cover, Black Wipeable, Size 3, Medium (8912432-82)	Back Cover, Black Wipeable, Size 3, Large (8912433-82)
Size 4	\$301	Back Cover, Black Wipeable, Size 4, Small (8912441-82)	Back Cover, Black Wipeable, Size 4, Medium (8912442-82)	Back Cover, Black Wipeable, Size 4, Large (8912443-82)

Replacement Black Wipeable Seat Cover

Additional Seat Cover only. Only available with Split Cushion.

Size	Price	Small	Medium	Large
Size 1	\$128	Seat Cover, Black Wipeable, Size 1, Small (8911411-82)	Seat Cover, Black Wipeable, Size 1, Medium (8911412-82)	Seat Cover, Black Wipeable, Size 1, Large (8911413-82)
Size 2	\$128	Seat Cover, Black Wipeable, Size 2, Small (8911421-82)	Seat Cover, Black Wipeable, Size 2, Medium (8911422-82)	Seat Cover, Black Wipeable, Size 2, Large (8911423-82)
Size 3	\$128	Seat Cover, Black Wipeable, Size 3, Small (8911431-82)	Seat Cover, Black Wipeable, Size 3, Medium (8911432-82)	Seat Cover, Black Wipeable, Size 3, Large (8911433-82)
Size 4	\$301	Seat Cover, Black Wipeable, Size 4, Small (8911441-82)	Seat Cover, Black Wipeable, Size 4, Medium (8911442-82)	Seat Cover, Black Wipeable, Size 4, Large (8911443-82)

Knee and Postural Supports

Seat Accessories



Anti-Thrust Wedges

Anti-Thrust Wedges (Pr.), Size 1 and 2 (8910400-2) (K0108)	\$56
Anti-Thrust Wedges (Pr.), Size 2 and 3 (8910400-4)(K0108)	\$56
Anti-Thrust Wedges (Pr.), Size 4 (8910400-12) (K0108)	\$77

Medial Thigh System



Medial Thigh System

Medial Thigh System, D: 3.5" x H: 3.5" x W: 3.5" (E0957,E1028) (8910408-1)	\$410
Medial Thigh System, D: 4.75" x H: 3.5" x W: 3.5" (E0957,E1028) (8910408-2)	\$410
Medial Thigh System, D: 6" x H: 3.5" x W: 3.5" (E0957,E1028) (8910408-3)	\$410

Abduction Supports Left and Right

Abduction Support w/ Cover

Abduction Supports Left and Right w/ Covers, Size 1 and 2 (8910402-10)	\$242
Abduction Supports Left and Right w/ Covers, Size 3 (8910402-20)	\$242
Abduction Supports Left and Right w/ Covers, Size 4 (8910402-40)	\$242

Hip Supports Left and Right w/ Covers

Hip Supports (Pair)

Hip Supports Left and Right w/ Covers, Size 1 and 2 (E0956) (8910404-10)	\$242
Hip Supports Left and Right w/ Covers, Size 3 (E0956) (8910404-20)	\$242
Hip Supports Left and Right w/ Covers, Size 4 (E0956) (8910404-40)	\$242

Hip Supports with PU Pad

Armrest, Single Post (8910619) is required

Hip Support Left and Right w/ PU Pad, 3.25" x 7.5" (E0956) (8910404-45)	\$657
---	-------



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Arm Supports

Armrest



Armrest with Cover

Armrest with Cover, Size 1, (Pair) (E0973) (8910615-2)	\$350
Armrest with Cover, Size 2, (Pair) (E0973) (8910616-2)	\$350
Armrest with Cover, Size 3, (Pair) (E0973) (8910617-2)	\$350
Armrest, Single Post, Size 4 (Pair) (E0973) (8910619)	\$350



Armrest Kits

Used for Armrests 1-3

Daily Adjustment Kit (8910615-11)	\$34
-----------------------------------	------

Upper Arm Support

Requires Fixed Side Supports



Upper Arm Support

Fixed Upper Arm Support with Cover, 5" x 4" (8910407-11)(K0108)	\$172
---	-------

Lateral Support Cushion and Fittings

Fixed Laterals

Fixed Side Supports w/ Covers (Pair)

Fixed Side Supports Left and Right w/ Covers, Type 1 (8910301-1)	\$242
Fixed Side Supports Left and Right w/ Covers, Type 2 (8910301-2)	\$242
Fixed Side Supports Left and Right w/ Covers, Type 3 (8910301-3)	\$242
Fixed Side Supports Left and Right w/ Covers, Type 4 (8910301-4)	\$242
Fixed Side Supports Left and Right w/ Covers, Type 5 (8910301-5)	\$242
Fixed Side Supports Left and Right w/ Covers, Type 6 (8910301-6)	\$242

Bars

Required selection for Swing Away Laterals. Quantity 2 required to mount left and right sides.

Straight Bars

Bar, Straight (81932) (Each)	N/C
------------------------------	-----

Offset Bars

Bar, Offset, 1" (81968) (Each)	N/C
Bar, Offset, 1.75" (82000) (Each)	N/C
Bar, Offset, 2.5" (81969) (Each)	N/C

Quantity:

Swing Away Lateral Support

Requires Bracket, Bar and Cushion selections

Swing-Away Lateral Support Bracket, Left (Each) (83900-L) (E1034)	\$267
Swing-Away Lateral Support Bracket, Right (Each) (83900-R) (E1034)	\$267

Lateral Cushions

Required selection for Swing Away Lateral Support. Sold in each, select quality 2 if wanting both left and right side.

Curved Lateral Cushion

Curved Cushion (Each), 3.5" x 7.5" (6175406) (E0956,E1034)	\$128
Curved Cushion (Each), 2.5" x 6.75" (6175405) (E0956,E1034)	\$128
Curved Cushion (Each), 2.5" x 6" (6175404) (E0956,E1034)	\$128

Standard Lateral Cushion

Standard Cushion (Each), 4" x 4" (6175402) (E0956,E1034)	\$128
Standard Cushion (Each), 2.2" x 4.5" (6175401) (E0956,E1034)	\$128

Protective Cushion

For lateral bar

Protective Cushion for Laterals, Black (89150) (Each)	\$112
---	-------

Quantity:



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Positioning Accessories

Positioning Belts

S-P: Shoulder to Pelvis and H-H: Hip to Hip

Torso Vest



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$232
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$232
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$232
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$232

Cross Vest



Cross Vest X-Small - S-P: 7.75" - 12.5" (Each) (85720-XS)(E0960)	\$220
Cross Vest Small - S-P: 10.25" - 13.5" (Each) (85720-S) (E0960)	\$220
Cross Vest Medium - S-P: 13.5" - 17.75" (Each) (85720-M) (E0960)	\$220
Cross Vest Large - S-P: 17.75" - 20.75" (Each) (85720-L) (E0960)	\$220
Cross Vest X-Large - S-P: 20.75" - 24.5" (Each) (85720-XL) (E0960)	\$220

4 Point Padded Hip Belt




4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$156
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$156
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M) (E0960)	\$156
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$156

2 Point Pelvic Belt



2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$98
---	------

Pelvic Harness



Pelvic Harness, Small - A:6.5" x B:7.5" (89108-0) (K0108)	\$284
Pelvic Harness, Medium - A:8.5" x B:8.75" (89108-1) (K0108)	\$284
Pelvic Harness, Large - A:10.5" x B:9.75" (89108-2) (K0108)	\$284

Headrest and Mounting Hardware

Mounting hardware required

Headrest	Price	Size
Large Laterals Headrest	\$332	Adjustable Width Headrest, Large Laterals, Black Cover, Short Mounting Hardware (9973792-85)
	\$250	Adjustable Width Headrest, Large Laterals, Black Cover, Long Mounting Hardware (9973793-85)
Large Laterals Headrest with Speakers	\$440	Adjustable Width Headrest, Large Laterals, with Speakers, Black Cover, Short Mounting Hardware (9971792-85)
		Adjustable Width Headrest, Large Laterals, with Speakers, Black Cover, Long Mounting Hardware (9971793-85)
Four Laterals Headrest	\$411	Adjustable Width Headrest, Four Laterals, Short Mounting Hardware (9975792-85)
		Adjustable Width Headrest, Four Laterals, Long Mounting Hardware (9975793-85)
Two Laterals Headrest	\$411	Adjustable Width Headrest, Two Laterals, Short Mounting Hardware (9977792-85)
		Adjustable Width Headrest, Two Laterals, Long Mounting Hardware (9977793-85)
Anatomic Head Support	\$271	Anatomic Head Support, Small (9923022-80) (E0955)
	\$281	Anatomic Head Support, Medium (9945022-80) (E0955)
	\$291	Anatomic Head Support, Large (9967022-80) (E0955)

Headrest Mounting System

Required selection. Required to mount headrest.

Headrest Mounting Hardware Multi-Adjustable, Short (98619-25) (E1033)	\$250
Headrest Mounting Hardware Multi-Adjustable, Long (98619-26) (E1033)	\$256



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Headrest Accessories

Headband System

Headband System for Anatomic Head Support Only (99500) (E0955) \$273

Headrest Mount

Used with other manufacturers' headrests

Universal Headrest Mount (87590) (E1028) \$105

Accessories

Trays

Tray requires mounting hardware

Two Rod Tray System

Tray, Two Rod System, Grey, Size 1 and 2, (Each) (8910725-1) (E0950) \$139
Tray, Two Rod System, Grey, Size 2 and 3, (Each) (8910725-2) (E0950) \$139
Tray, Two Rod System, Grey, Size 4, (Each) (8910725-3) (E0950) \$139

Two Rod Tray System Hardware

Swing Away Removable Hardware, Size 1 and 2 (Each) (X8910725-12) (E1028) \$230
Swing Away Removable Hardware, Size 2 and 3 (Each) (X8910725-22) (E1028) \$230
Swing Away Removable Hardware, Size 4 (Each) (X8910725-32) (E1028) \$230

Tray Accessories

Not available with Size 4

Toy Bar for Tray (8910740) \$272
Arm Pads, Mounts on Top of Tray, Size 1, 2, 3 (8910730) \$350

Canopy and Hardware

Canopy

Canopy, Black, Size 1 and 2 (8740110) \$398
Canopy, Black, Size 3 and 4 (8740120) \$398

Canopy Hardware

Holders for Canopy (9589100) \$56

Footrests for Bases

Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Foot Supports for Bases

Select footrest type and desired mounting bar

Footrest - Angle and Height Adjustable

Angle and Height Adjustable Foot Plate, Small, 12" x 7.5" (95145-1) (K0040) \$325
Angle and Height Adjustable Foot Plate, Medium, 14" x 9" (95145-2) (K0040) \$325
Angle and Height Adjustable Foot Plate, Large, 16" x 9.5" (95145-3) (K0040) \$325

Footrest - Multi-Adjustable

Mounting hardware included. Do not quote mounting bar. Not available with a Foot Box.

Multi-Adjust Foot Support, Small - 6" - 10.5" (8910852-1) \$495
Multi-Adjust Foot Support, Medium - 6" - 14.5" (8910852-2) \$495
Multi-Adjust Foot Support, Large - 6" - 17.5" (8910852-3) \$495

Footrest Mounting Bar

Required for Angle and Height Adjustable foot plate

Short Mounting Bar - 6-10.5" (8910821-1) N/C
Medium Mounting Bar - 6-14.5" (8910821-2) N/C
Long Mounting Bar - 6-17.5" (8910821-3) N/C

Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjustable (8910855) \$82



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Foot Box Kits

Requires Angle and Height Adjustable Foot Plate above and rod. Select desired size Foot Box and corresponding Foot Plate padding to match footrest height. Not available with Multi-Adjustable Footrests.

Foot Box Kit

Foot Box Kit, Small Plate (95145-1BOX)	\$240
Foot Box Kit, Medium Plate (95145-2BOX)	\$255
Foot Box Kit, Large Plate (95145-3BOX)	\$275

Strap for Foot Box

Strap for Foot Box (95145-20)	\$98
-------------------------------	------

Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD)	\$120
Foot Plate Padding, Size 2 (95145-2PAD)	\$165
Foot Plate Padding, Size 3 (95145-3PAD)	\$205

Padding for Foot Plate Rod

Padding for Foot Plate Rod (95145-150)	\$64
--	------

Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$165
Neoprene Sandals - Medium (85696-2)	\$165
Formed Shoeholder - Small (86147)	\$195
Formed Shoeholder - Medium (86148)	\$195
Footstraps w/ fix lock (85096) (E0951)	\$42
Ankle Straps, X-Small - C: 6.25"-9" (85745-XS)	\$202
Ankle Straps, Medium - C: 9"-12.25" (85745-M)	\$202
Ankle Straps, X-Large - C: 11.5"-15" (85745-XL)	\$202



Attaches to mounting bar

Calf Support for Centrally Mounted Foot Plate (85131)	\$245
---	-------

High-Low:xo Indoor Base

High-Low:xo Indoor Bases

Manual Lifts

High-Low:xo Manual Lift, Indoor, Size 1 (89010522) (K0108)	\$2196
High-Low:xo Manual Lift, Indoor, Size 2 (89020522) (K0108)	\$2356
High-Low:xo Hydraulic Lift, Indoor, Size 3 (89030542) (K0108)	\$3537

Power Lifts

High-Low:xo Indoor, Power Lift, Size 3 (89030532) (K0108)	\$4124
---	--------

High-Low:xo Indoor Bases Accessories

Adapter Mounting Plate

Required for size 4 x:panda on a High-Low:xo, Size 3

Adapter Plate for x:panda sz 4 on High-Low:xo 3 (8911610-99)	\$211
--	-------

Push Brace

Required selection for indoor High-Low:xo Base

Push Brace, Fixed Angle (890009-04)	\$194
Push Brace, Angle Adjustable (8910555)	\$465

Anti-Tip

Anti-Tip for Indoor High-Low:xo Indoor (890007)	\$248
---	-------

Chilla Mobility Base

Chilla Frame

Chilla Folding Frame

Only available with size 1 and size 2

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234	\$2630
--	--------

Shape Adapter

x:panda Shape Adapter Plate (881091) K0208	\$875
--	-------

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089)	\$224	Vent Tray (881095)	\$330
Under Storage Basket (881092)	N/C	O2 Holder (881096)	\$175
IV Pole (881097)	\$475		



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

x:panda

Multi Frame Mobility Base

R82 Multi Frame

Multi Frame Configuration

Multi Frame Size 3, Silver, 12.5" Solid Tire, Std. Wheel Locks (973140210)	\$2995
Multi Frame Size 3, Silver, 12.5" Solid Tire, Drum Brake (973141210)	\$2995
Drum Brake for Multi Frame	\$600
<i>Required selection with Drum Brake, Multi Frame configuration</i>	

Multi Frame Accessories

Adapter

Adapter for x:panda on Multi Frame, Size 1-3 (8911610-98)	\$210
Adapter for x:panda on Multi Frame, Size 4 (8911610-99)	\$211

Accessories

Anti-Tips (970957) (E0971)	\$165
Transit (970905)	\$240
Vent Tray for Multi Frame (970815) (E1029)	\$525

Ki Mobility Bases

Separate Order form for base must be filled out. Order form is located at Kimobility.com. Footrests are not provided with the base. If footrest are desired, selections must be made above in the "Mobility Base Footrest" section.

Ki Mobility Bases

Tyke x:panda

Only available with Size 1 and Size 2. Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Tyke for x:panda	N/C
------------------	-----

Little Wave Flip for x:panda

Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Flip for x:panda	N/C
------------------	-----