



by Etac
Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Infant Seat

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

		Infant	1	2
	Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
	Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
	Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
	Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
	Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
	Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
	Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
	Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
	Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Seat Bases

Please see measurement guide for growth information. A selection must be made in each section below.

Set Base

Includes base and plate

Size 1 Infant Base

See chart above for size details

\$895

Recline Options

Selection Required

Adjustable Stepless Recline, Toolless (K0108)

\$443

Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (K0108)

\$723

Back

Upper Back with Quick Release Mount

\$260

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0625.5 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Infant Seat

Cushions and Covers

Seat Cushion

Cushion Color

Cushion including seat, back, seat pit, upper arm cushions. Order covers for side support belt separately. Required Selection.

Pink	N/C
Light Blue	N/C
Black Wipeable	\$21

Extra Cushion and Cover

An extra cushion and cover set, in addition to the installed cushion and cover selected above. Comes uninstalled.

Pink	\$339
Light Blue	\$339
Black Wipeable	\$364


Extra Cover

Select for an extra cover (no cushion). Comes uninstalled. Optional selection.

Pink	\$183
Light Blue	\$183

Headrests and Hardware Options

Headrest Style



Type of Head Support

Includes Short Rod length. Includes Cushion.

Curved, Small- Black Wipeable - W: 7", H: 5" (E0955)	\$220
Curved, Small -Black - W: 7", H: 5" (E0955)	\$220
Anatomic, Small A: 5.5" x B: 5.5" x C: 9.5" - Black, Wipeable (E0955)	\$271
Anatomic, Small A: 5.5" x B: 5.5" x C: 9.5" - Black (E0955)	\$271

Headrest Mounting Hardware

2 pivot points, standard 80mm rod, height adjustment rod 905037

Standard Mounting Hardware (E1033)	\$245
------------------------------------	-------



Headrest Options

Knob replaces bolt for tool-less adjustment of headrest angle

Daily Adjustment Kit Standard - Set of Two	\$32
Daily Adjustment Kit Standard - Set of Three	\$47

Vest, Straps and Belts



Side Support Belt

Infant Side Support Belt, Black Wipeable (E0978)	\$219
Infant Side Support Belt, Black (E0978)	\$219

Side Support Extra Cover

Extra Side Support Cover, Black, Wipeable	\$99
Extra Side Support Cover, Black	\$99

Cross Vest

Cross Vest X-Small - S-P: 7.75" - 12.5" (E0960)	\$220
Cross Vest Small - S-P: 10.75" - 13.5" (E0960)	\$220

Torso Vest

Torso Vest X-Small - S-P: 7.75" - 10.75" (E0960)	\$232
Torso Vest Small - S-P: 10.75" - 12.5" (E0960)	\$232

Padded Hip Belt

4-Point Padded Hip Belt, 1" X-Small - H-H: 6.75" - 12.5"(E0978)	\$156
4-Point Padded Hip Belt, 1" Small - H-H: 8.25" - 13.75" (E0978)	\$156
2-Point Pelvic Belt Black - L: 13-14.5" (E0978)	\$98



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Trays and Accessories

Canopy

Black Canopy - Size 1 _____ \$339



Tray

Soft Tray for Infant Seating, Black - W: 12.5" x L: 7.5" _____ \$156

Canopy Brackets

Canopy Brackets with Push Handle for Size 1 _____ \$203

Footrests for Mobility Bases

Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrests

Footrest Type

Angle and Height Adjustable Foot Plate, Small, 12" x 7.5" (95145-1)(K0040) _____ \$325

Mounting Bar

Required selection with Angle and Height Adjustable Foot Plate

Short Mounting Bar - 6-10.5" (8910821-1) _____ N/C

Footrest Accessories

Neoprene Sandals - Small _____ \$165
Neoprene Sandals - Medium _____ \$165

Footstraps w/ fix lock (E0951) _____ \$43
Ankle Straps, X-Small - C: 6.25" - 9" _____ \$197

Chilla Mobility Base

Shape Adapter required

Chilla Frame

Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234 _____ \$2630

Shape Adapter

x:panda Shape Adapter Plate (881091) K0208 _____ \$875

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) _____ \$224

Vent Tray (881095) _____ \$330

IV Pole (881097) _____ \$475

O2 Holder (881096) _____ \$175

High-Low:xo Base and Options

High-Low:xo Bases and Options

Base

Manual High-Low:xo bases require a Push Brace selection

Manual Lift - Indoor Size 1 (K0108) _____ \$1976

Indoor- Power Lift, Size 1 (K0108) _____ \$3175

Anti-Tips for High-Low

Anti-Tips - Indoor Base (890007-1) _____ \$223

Push Brace for High-low:xo - Indoor

Required with High-Low:xo

Fixed Angle Push Brace _____ \$198

Angle Adjustable Push Brace _____ \$419

Ki Mobility Bases

Separate order form for base must be filled out. Order form is located at www.kimobility.com. Footrests are not provided with base. If footrests are desired, please make selection at Footrests for Mobility Bases section above.

Tyke for x:panda _____ N/C

Tyke size 1 is only available for Infant Package, Tyke size 2 is only available for Shape Simple/Advanced size 1

Flip for x:panda _____ N/C

Shape size 2 only for Flip