



by Etac
Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

x:panda Shape - Simple Seat

	Infant	1	2
 Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
 Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
 Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
 Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
 Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
 Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
 Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Shape Base Options

Please see measurement guide for growth information. A selection must be made in each section below.

Shape Base Options

Selection required

Seat Base

Size 1 Simple Seat (E2292)	\$931
Size 2 Simple Seat (E2292)	\$931

Seat Plate

Selection required

Seat Plate Size 1

Size 1, Seat W: 10", Lower Back W: 8.75"	N/C
Size 1, Seat W: 11.75", Lower Back W: 11.75"	N/C

Seat Plate Size 2

Size 2, Seat W: 12.5", Lower Back W: 12.5"	N/C
Size 2, Seat W: 15.75", Lower Back W: 15.75"	N/C

Recline Options

Selection required

Back Recline

Adjustable Stepless Recline, Toolless (K0108)	\$443
Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (K0108)	\$723

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0627.5 Pricing and specifications subject to change without notice.



by Etac
Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Back and Cushions

Upper Back Size

Select upper back base then select Upper and Lower back cushion and cover. Includes quick release detachment (#9050374)



Size 1 Base

Size 1 Base Upper - H: 7" (905073) (E2291) \$260
Height range: 12.25" - 16.25"
Size 1 Base Upper - H: 8.75" (905074) (E2291) \$260
Height range: 13.75" - 18.50"

Size 2 Base

Size 2 Base Upper - H: 7" (905073) (E2291) \$260
Height range: 13.5" - 17.25"
Size 2 Base Upper - H: 8.75" (905074) (E2291) \$260
Height range: 15" - 19.25"
Size 2 Base Upper - H: 10.75" (905075) (E2291) \$260
Height range: 17" - 21.75"

Upper Back Cushion and Cover

Upper back cushion height must match Upper Back Base height

Upper Back Cushion Size

Upper Back Cushion-W: 7"xH: 5.5" (902101) N/C
Upper Back Cushion-W: 7"xH: 7" (902103) N/C
Upper Back Cushion-W: 8.75"xH: 7" (902104) N/C
Upper Back Cushion-W: 10.25"xH: 7" (902105) N/C
Upper Back Cushion-W: 8.75"xH: 8.75" (902106) N/C
Upper Back Cushion-W: 10.25"xH: 8.75" (902107) N/C
Upper Back Cushion-W: 11.75"xH: 8.75" (902108) N/C
Upper Back Cushion-W: 13.25"xH: 8.75" (902109) N/C
Upper Back Cushion-W: 10.25"xH: 10.75" (902110) N/C
Upper Back Cushion-W: 11.75"xH: 10.75" (902111) N/C
Upper Back Cushion-W: 13.75"xH: 10.75" (902112) N/C
Upper Back Cushion-W: 15"xH: 10.75" (902113) N/C

Upper Back Cushion Cover

Selected required

Red (-50) N/C
Grey (-60) N/C
Black Wipeable (-70) \$21

Lower Back Cushion and Cover

The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.

Lower Back Cushion Size

Shape Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120) N/C
Shape Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902121) N/C
Shape Size 1 Lower Back Cushion - W: 11" x H: 5.5" (902122) N/C
Shape Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123) N/C
Shape Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124) N/C
Shape Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125) N/C
Shape Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126) N/C

Lower Back Cushion Cover

Selection required

Red (-50) N/C
Grey (-60) N/C
Black Wipeable (-70) \$21

Daily Adjustment Kit - Upper Back

Allows for toolless adjustable of back height. Ideal for frequent adjustment.



Daily Adjustment Kit

Daily Adjustment Kit - Upper Back (905081) \$21

Spacer Blocks

Standard Spacer Block Holder without Cushion

Cushion widths must match Spacer Block Holder width

Spacer Block Holder, Size 1 - Back widths 7", 8.75", 10.25" (903100-1) \$183
Spacer Block Holder, Size 2 - Back widths 8.75", 10.25", 11.75", 13.5" (903100-2) \$183

Standard Spacer Block Cushion Cover

Black (-80) N/C
Black Wipeable (-70) N/C

Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011) N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012) N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013) N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014) N/C
Spacer Block Cushion Size 5 - W: 13.25" x H: 2.25" (903015) N/C



Cushions and Covers

Seat Cushions

Seat Cushion Size - One Piece

One Piece Cushion, Size 1 - W: 9.75" x L: 19.75" (901101) N/C
One Piece Cushion, Size 1 - W: 11.75" x L: 19.75" (901103) N/C
One Piece Cushion, Size 2 - W: 12.50" x L: 25.5" (901105) N/C
One Piece Cushion, Size 2 - W: 15.75" x L: 25.5" (901107) N/C

Seat Cushion Size - Split

Required when Quick Release Pommel is selected

Split Cushion, Size 1 - W: 10.25" x L: 19.75" (901114) N/C
Split Cushion, Size 1 - W: 11.75" x L: 19.75" (901115) N/C
Split Cushion, Size 2 - W: 12.5" x L: 25.5" (901117) N/C
Split Cushion, Size 2 - W: 15.75" x L: 25.5" (901119) N/C

Seat Cover

Required selection for cover color

Red (-50) N/C
Grey (-60) N/C
Black Wipeable (-70) \$21

Headrest and Hardware Options

Headrest and Hardware Options

Includes Cushion. Select desired headrest in chart and choose the corresponding required hardware below.

Style	Size	Price	Type
Curved	Small	\$220	Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)
			Curved, Small - Black - W: 7", H: 5" (E0955)
	Medium	\$240	Curved, Medium - Black Wipeable W: 10.2", H: 6.3" (E0955)
			Curved, Medium - Black W: 10.2", H: 6.3" (E0955)
Anatomic	Small	\$271	Anatomic, Small - Black Wipeable A: 5.5" x B: 5.5" x C: 9.5" (E0955)
			Anatomic, Small - Black A: 5.5" x B: 5.5" x C: 9.5" (E0955)
	Medium	\$281	Anatomic, Medium - Black Wipeable A: 6" x B: 5.5" x C: 10" (E0955)
			Anatomic, Medium - Black A: 6" x B: 5.5" x C: 10" (E0955)
	Large	\$291	Anatomic, Large - Black Wipeable A: 6.5" x B: 5.5" C: 12.5" (E0955)
			Anatomic, Large - Black A: 6.5" x B: 5.5" C: 12.5 (E0955)
Adjustable		\$332	Adjustable w/ Large Laterals - Black A: 4.75" x B: 6 x C: max 9.5" (E0955)
		\$440	Adjustable with Large Laterals and Speakers - Black A: 4.75" x B: 6: max 9.5" (E0955)
		\$411	Adjustable with Four Laterals - Black A: 4.75" x B: Max 9.5" C: max 8" (E0955)



by Etac

Phone (877) 227-0097 | Fax (704) 882-0751

sales.us@etac.com

x:panda shape - Simple Seat

Mounting Hardware for Curved and Anatomic Headrest

Headrest Mounting Hardware

Mounting Brackets, Standard (E1033)	\$245
Mounting Brackets, Extended (E1033)	\$255
Mounting Brackets, 3D Adjustment (E1033)	\$266

Mounting Hardware for Adjustable Headrest Style

Short Mounting Hardware (E1033)	\$250
Long Mounting Hardware (E1033)	\$250

Headrest Options



Daily Adjustment Kits

Required with Curved and Anatomic Head support

Daily Adjustment Kit Standard - Set of Two (905034)	\$32
Daily Adjustment Kit Extended - Set of Three (905035)	\$47

Additional Headrest Options

Headband System	\$250
Universal Headrest Mount	\$67

Arm Supports

Arm Support

Arm Supports

Cover not included, select separately

Upper Arm Support, Fixed W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$339
Upper Arm Support, Fixed W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$339

Arm Supports Cover

Black (-80)	\$79
Black Wipeable (-70)	\$79

Daily Adjustment Kit

For toolless height adjustment of arm supports

Height Adjustment Kit of Arm Supports (8910615-11)	\$21
--	------

Upper Arm Supports - Fixed

Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.



Upper Arm Supports, Fixed, Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$260
<i>For size 1 Base</i>	
Upper Arm Supports, Fixed Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$260
<i>For Size 2 Base</i>	

Upper Arm Supports - Adjustable

Not compatible with Swing-Away Lateral Supports

Upper Arm Supports, Adjustable, Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$349
<i>For Size 1 Base</i>	
Upper Arm Supports, Adjustable, Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$349
<i>For size 2 Base</i>	

Knee and Postural Supports

Abduction and Adduction



Outer Knee Support

Outer Knee Support, Standard Set for Size 1, Black - L: 6", H: 3.25" (905005-70) (E0956)	\$219
Outer Knee Support, Standard Set for Size 1, Black Wipeable - L: 6", H: 3.25" (905005-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black - L: 7.75", H: 3.5" (905006-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black Wipeable - L: 7.75", H: 3.5" (905006-70) (E0956)	\$219

Quick Release Pommel

Quick Release Pommel For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$364
Quick Release Pommel For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905082-80) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905082-70) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905078-80) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905078-70) (E0957)	\$364

Postural Supports

Width Adjustable Pelvis Supports

Mounted from the back



Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" (905041-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" (905041-70) (E0956)	\$245
Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" (905042-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Large Black	\$245

Width Adjustable Hip and Pelvis Supports

Mounted from the sides



Width Adjustable Hip Supports, Small Black - L: 4" H: 4.25" (905044-80) (E0956)	\$349
Width Adjustable Hip Supports, Small Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$349
Width Adjustable Hip Supports, Large Black - L: 5.5" H: 5" (905045-80) (E0956)	\$349
Width Adjustable Hip Supports, Large Black	\$349

Trays and Accessories

Trays and Activity Centers

Tray



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950) \$334
For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.



Shape 6 Tray with Edge Size 2 (8910725-2) (E0950) \$334
For Size 2. Double rod, Grey ABS.



Shape 8 Tray with Edge Size 1 (8976350) (E0950) \$334



Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028) \$443

Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028) \$312

Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028) \$443

Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028) \$312

Additional Tray Accessories



Toy Bar (8910740) \$250
For shape 5 and 6

For Shape 5 and 6

Arm Pads for Tray (8910730) \$315

Canopy and Mounting Hardware

Canopy does not come with brackets. Both the canopy and bracket must be quoted when mounting to the shape.



Canopy Without Brackets

Black Canopy - Size 1 (8740110) \$339

Black Canopy - Size 2 (8740120) \$339



Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1) \$203

Canopy Brackets with Push Handle for Size 2 (905039-2) \$203

Accessory Holders and Brackets

Remote Control Holder Size 1 - Arm Support Mounted (905030-1) \$68

Remote Control Holder Size 2 - Arm Support Mounted (905030-2) \$68

Remote Control Holder - Back Mounted (905030) \$131

Bracket for Communication Aid - Arm Support Mounted (905033-1) \$82

Vests and Belts

More options can be found at etac.com

Vests

Torso Vests



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960) \$232

Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960) \$232

Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960) \$232

Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960) \$232

Cross Vests



Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960) \$220

Cross Vest Small - S-P 10.25" - 13.5" (85720-S) (E0960) \$220

Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960) \$220

Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960) \$220

Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960) \$220

Fixlocks Bracket



Higher Brackets (8910260) \$77
Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.

Wider Brackets (905068) \$57



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Pelvic Belts and Harnesses

4-Point Belts

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$156
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$156
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M) (E0960)	\$156
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$156

Pelvic Harness

Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0)	\$261
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1)	\$261
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2)	\$261

Pelvis Belt

Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$175
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-S) (E0978)	\$175
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$175

2-Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$98
---	------

Lateral Support Cushion and Fittings

if choosing Laterals, both the left and right cushions must be selected

Left Lateral Supports

Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Right Lateral Supports

Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Fitting for Lateral Supports

Selection required, fittings are sold as a pair

Type of Fitting

Fixed Fitting (E1028)	\$172
Swing-Away and Angle Adjustable Fitting (E1028)	\$260

Footrests for Bases

Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrest

Select footrest type and desired mounting bar. Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrest - Angle and Height Adjustable

Angle and Height Adjustable Foot Plate, Small, 12" x 7.5" (95145-1) (K0040)	\$325
Angle and Height Adjustable Foot Plate, Medium, 14" x 9" (95145-2) (K0040)	\$325
Angle and Height Adjustable Foot Plate, Large, 16" x 9.5" (95145-3) (K0040)	\$325

Footrest Mounting Bar

Required for Angle and Height Adjustable Foot Plates

Short Footrest Mounting Bar - 6-10.5"	N/C
Medium Footrest Mounting Bar - 6-14.5"	N/C
Long Footrest Mounting Bar - 6-17.5"	N/C

Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjustable (8910855)	\$82
--	------

Individual Multi-Adjustable Foot Supports

Mounting hardware comes standard. Not compatible with the Footbox.

Individual Multi-Adjust Foot Supports, Small - 6" - 10.5" (8910852-1)	\$495
Individual Multi-Adjust Foot Supports, Medium - 6" - 14.5" (8910852-2)	\$495
Individual Multi-Adjust Foot Supports, Large - 6" - 17.5" (8910852-3)	\$495



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Foot Box Kits

Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi-Adjustable Footrest.

Foot Box Kit

Required selection

Foot Box Kit, Size 1 Plate (95145-1BOX)	\$240
Foot Box Kit, Size 2 Plate (95145-2BOX)	\$255
Foot Box Kit, Size 3 Plate (95145-3BOX)	\$275

Strap for Foot Box

Required with all sizes

Strap for Foot Box (95145-20)	\$98
-------------------------------	------

Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD)	\$120
Foot Plate Padding, Size 2 (95145-2PAD)	\$165
Foot Plate Padding, Size 3 (95145-3PAD)	\$205

Padding for Foot Plate Rod

Required selection

Padding for Foot Plate Rod (95145-150)	\$64
--	------

Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$165
Neoprene Sandals - Medium (85696-2)	\$165
Formed Shoeholder - Small (86147)	\$195
Formed Shoeholder - Medium (86148)	\$195
Footstraps w/ fix lock (85096) (E0951)	\$43
Ankle Straps, X-Small - C: 6.25" - 9" (85745-XS)	\$197
Ankle Straps, Medium - C: 9" - 12.25" (85745-M)	\$197
Ankle Straps, X-Large - C: 11.5" - 15" (85745-XL)	\$197



Calf Support for Centrally Mounted Foot Plate (85131)	\$245
Attaches to Mounting Bar	

Chilla Base

Chilla Frame

Adaptor for x:panda Shape required

Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234	\$2630
--	--------

Shape Adapter

Required to mount to Chilla Base

x:panda Shape Adapter Plate (881091) K0208	\$875
--	-------

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089)	\$224
Under Storage Basket (881092)	N/C
Vent Tray (881095)	\$330
O2 Holder (881096)	\$175
IV Pole (881097)	\$475

High-Low:xo Base and Options

High-Low:xo Base

Manual Lifts

Manual High-Low bases require a Push Brace selection

High-low:xo Frame Manual Lift, Indoor, Size 1	\$1976
High-low:xo Frame Manual Lift, Indoor, Size 2	\$2121

Outdoor High-Low:xo

High-low:xo Frame Outdoor, Size 2, Manual Lift	\$3329
High-low:xo Frame Outdoor, Size 2, Power Lift	\$4266

Power Lifts

High-low:xo Frame Power Lift, Indoor, Size 1	\$3175
High-low:xo Frame Power Lift, Indoor, Size 2	\$3409

High-Low:xo Options

Transit - Outdoor High-Low:xo

Transit (890005) (E1022)	\$260
--------------------------	-------

Push Brace

Available only with Manual Indoor Lifts

Fixed Angle Push Brace (890009-04)	\$198
Angle Adjustable Push Brace (8910555)	\$419

Anti-Tips on the High-Low:xo

Anti-Tips - Indoor Base (890007-1)	\$223
Anti-Tips - Outdoor Base (890008)	\$240



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Multi Frame Mobility Base and Options

Adaptor plate for x:panda Shape required

Multi Frame		
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel		
Locks	\$2995	
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake	\$3595	
Adaptor for x:panda on Multi Frame - for x:panda Shape 1 & 2	\$210	
Base options		
Anti-tips		\$165
Transit for Multi Frame		\$240
Vent Tray		\$525
O2 Holder		\$395

Ki Mobility Bases

Separate order form for base must be filled out. Order form located at Kimobility.com. Footrests are not provide with base. If footrests are desired, make a selection in the Footrest for Base, section.

Ki Mobility Bases		
Price for base calculated with Ki Mobility's order form. please refer to Tyke order form to verify compatibility.		
Flip for x:panda		N/C
Shape size 2 only for Flip		
Tyke for x:panda		N/C
Tyke size 1 is only available for Infant Package, Tyke size 2 is only available for Shape Simple/Advanced size 1		