



MYCO
MEDICAL

TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET

REDI-CUT®


Disposable Safety Scalpel with TECHNOCUT® Premium Blade and Retractable Sheath

Portability and surgical quality are joined together in one convenient product with the REDI-CUT® Disposable Safety Scalpel, a disposable scalpel that boasts a TECHNOCUT® Premium stainless steel blade and an easy-to-use retractable safety sheath. Available sterile or non-sterile.

- Top activation allows for ambidextrous use
- Audible “clicks” confirm sheath positions – “Safe,” “Ready to use,” and “Final lock”
- Scalpel blades produced from high-quality Swedish stainless steel (TECHNOCUT® Premium blade)

Pictured: REDI-CUT® Disposable
Safety Scalpel, #10 Blade

- Sheath retracts fully for clear line of sight and access to surgical blade
- Blade can be covered for safe passing
- Final lock for safe handling and disposal
- Features non-slip, ergonomic handle with integrated metric ruler and depth measure
- Scalpels manufactured in an ISO 13485 certified facility



Sheath retracts fully for
clear line of sight



Final lock for safe handling
and disposal

BRAND (Company)	Top Activated	Final Lock	Audible Clicks
REDI-CUT® Disposable Safety Scalpel with Retractable Sheath	✓	✓	✓
Bard-Parker™ Disposable Safety Scalpel	✓		✓
Shackelford (Deroyal)		✓	✓
Futura (Merit Medical)			✓
Safety Scalpel (Medline)	✓		✓
Safety Scalpel (Carefusion)		✓	✓
Southmedic Plus (Southmedic)	✓	✓	✓
TECHNOCUT® Disposable Safety Scalpel (MYCO Medical)	✓	✓	✓



Sterile REDI-CUT® Disposable Safety Scalpels with Retractable Sheath

BLADE SIZE	MYCO CODE	UNIT OF MEASURE
#10	RCRSS10	10/bx, 100bx/cs
#11	RCRSS11	10/bx, 100bx/cs
#15	RCRSS15	10/bx, 100bx/cs

Non-Sterile REDI-CUT® Disposable Safety Scalpels with Retractable Sheath are available upon request. Minimums and lead times apply.

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com

www.mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2020 MYCO Medical, Inc. All rights reserved.

LIT-00660 Rev 0 Eff 01.08.2020