



INTRODUCING THE **NEW** **PEDIATRIC FLOWMETER** FROM RHYTHM HEALTHCARE

Rhythm Healthcare's **NEW** Pediatric Flowmeter is a convenient, lightweight and easy to use device. Using a standard 5L Concentrator, you can meet the needs of your low flow pediatric patients in the home setting.

- ✓ Delivers flow rates as low as 1/16LPM
- ✓ Just 1.50 lbs.
- ✓ Requires no electricity
- ✓ Designed to run from a Standard Stationary Oxygen Concentrator such as the Rhythm Healthcare LM5BA
- ✓ Simple Operation
- ✓ Easy to Read flowmeter



RHYTHM
HEALTHCARE

Item # **FMP01**

PEDIATRIC FLOWMETER INFORMATION

PRODUCT SPECIFICATIONS

Input pressure rating	34.5 ± 3.45 kPa (5 ± 0.5 psi)
Oxygen input:	2 L/min
Oxygen output:	0 - 750 cc/min (0 - ¾ L/min)
Flow adjustable in increments	50 cc/min (1/16 L/min)
Flow accuracy	± 10% of full scale
Size (L x W x H):	4.72 x 4.76 x 5.93 in (120 x 121 x 150.5 mm)
Net Weight:	1.50 lbs (0.68 kg)
Shipping weight:	2.54 lbs (1.15 kg)
Operating Temperature:	50° - 104° F (10° - 40° C)
Storage & Transportation Temperature:	-40° - 158° F (-40 - 70° C)



INCLUDED WITH PEDIATRIC FLOWMETER

Description	Quantity
Pediatric Flowmeter with Integrated Stand	1
Crushproof Tubing 25ft	1
Crushproof Tubing 7ft	2
Connector	1
Humidifier Bottle Adapter	1
Instructions	1

