

KNEE WALKER W/ SHOCK ABSORBERS & TOTE BAG

300 lb Weight Capacity



ITEM #	DESCRIPTION
KN2100-BK	Knee Walker (Black)
KN2100-OR	Knee Walker (Orange)

FEATURES

- 10" Front Wheels & 8" Rear Wheels
- Detachable Fabric Tote Bag
- Height Adjustable Handlebars
- Knee Pad is Height & Depth Adjustable
- Shock Absorbers for Smooth Movement
- Multi-functional Brakes
- Folds into Compact Size for Easy Storage/Transportation

SPECIFICATIONS

Product Weight:	28.5 lbs
Frame:	Steel
Knee Pad Length:	12.7"
Handle Height:	33.25" - 40"
Handlebar Length:	16.9"
Overall Dimensions:	32.4" - 35.4"L x 17.5" W x 33.25" H
Weight Capacity:	300 lbs

WARRANTY

Rhythm Healthcare, LLC offers a limited 1 year warranty from defects in materials and craftsmanship of the product to the original purchaser. This warranty does not cover normal wear and tear or misuse of the product. Rhythm Healthcare, LLC reserves the right to determine if a warranted product should be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare, LLC dealer.

V2: 8/16/24