

# P2 PORTABLE OXYGEN CONCENTRATOR BATTERY



| ITEM #    | DESCRIPTION   |
|-----------|---|
| P2BY001-2 | P2/P2-E6/P2-E7 Portable Oxygen Concentrator Lithium-ion Battery |

- Compatible with P2/P2-E6/P2-E7

## SPECIFICATIONS

|                         |   |
|-------------------------|---|
| Rechargeable Battery    | 14.4V   |
| Battery Rated Capacity: | 6.8Ah   |
| Battery                 | 6800mAh, duration up to 5 hours on setting 1.<br>Approx. 4 hours charge time with AC. |

## WARRANTY

1 Year