





|                                 |                                |
|---------------------------------|--------------------------------|
| <b>Commercial product name:</b> | <b>SED-CHEK®2 ESR CONTROLS</b> |
|---------------------------------|--------------------------------|

**4. FIRST AID MEASURES:**

General information: no special measures required  
After inhalation supply fresh air, consult doctor in case of complaints after skin contact  
Immediately wash with water and soap and rinse thoroughly  
Generally the product does not irritate the skin.  
After eye contact, rinse opened eye for several minutes under running water. Then consult a doctor.  
After swallowing rinse mouth with water. Seek medical attention and appropriate follow-up information for doctor.  
The following symptoms may occur: Skin and eye irritation.  
Medical conditions aggravated by exposure: None expected.

**5. FIRE-FIGHTING MEASURES:**

|                      |   |
|----------------------|---|
| Extinguishing media: | CO <sub>2</sub> extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant form. |
|----------------------|---|

**6. ACCIDENTAL RELEASE MEASURES:**

|                 |   |
|-----------------|---|
| After spillage: | Measures for cleaning/collecting:<br>Absorb liquid components with liquid-binding material.<br>Pick up mechanically.<br>Clean the affected area carefully. Suitable cleaners are disinfectant.<br>Dispose contaminated material as waste.<br>Additional information: no dangerous substances are released |
|-----------------|---|

**7. HANDLING AND STORAGE:**

|                        |   |
|------------------------|---|
| 7.1 Handling:          | Information for safe handling: no special precautions are necessary if used correctly.<br>Information about protection against explosions and fires: No special measures required |
| 7.2 Storage            | Further information about storage conditions: none<br>Requirements to be met by storerooms and receptacles: no special requirements   |
| Do not store together: | No restrictions.  |



|                                 |                                |
|---------------------------------|--------------------------------|
| <b>Commercial product name:</b> | <b>SED-CHEK®2 ESR CONTROLS</b> |
|---------------------------------|--------------------------------|

**8. EXPOSURE CONTROLS/PERSONAL PROTECTION:**

|                         |   |
|-------------------------|---|
| TLV:                    | Threshold Limit Value established by ACGIH: substance not listed. |
| Respiratory protection: | Not required  |
| Eye protection:         | Safety glasses  |
| Hand protection:        | Protective gloves   |

**9. PHYSICAL AND CHEMICAL PROPERTIES:**

|   |                 |
|---|-----------------|
| Physical state:                                   | LIQUID REAGENTS |
| Color:  | Red             |
| Odor:   | Light           |
| pH value at 20°C:                                 | 7.0-9.0         |
| Boiling point                                     | N/A             |
| Melting point:                                    | N/A             |
| Flash point:                                      | N/A             |
| Substance does not have any oxidising properties: | N/A             |
| Ignition temperature                              | N/A             |
| Explosion limits:                                 | N/A             |
| Vapor pressure:                                   | N/A             |
| Density:  | N/A             |
| Solubility in water:                              | Fully           |
| Viscosity   | N/A             |

**10. STABILITY AND REACTIVITY:**

|                                   |                               |
|-----------------------------------|-------------------------------|
| Hazardous Reactions:              | None when used appropriately. |
| Hazardous decomposition products: | N/A.                          |



|                                 |                                |
|---------------------------------|--------------------------------|
| <b>Commercial product name:</b> | <b>SED-CHEK®2 ESR CONTROLS</b> |
|---------------------------------|--------------------------------|

**11. TOXICOLOGICAL INFORMATION:**

|                                       |   |
|---------------------------------------|---|
| Acute toxicity:                       | Primary irritant effect: On the skin: No irritant effect  |
| On the eye:                           | No irritant effect sensitization: No sensitizing effect known   |
| Additional toxicological information: | The product is not subject to classification according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version.<br><br>When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us |
| Target organs:                        | Not applicable  |

**12. ECOLOGICAL INFORMATION:**

|                     |  |
|---------------------|--|
| Water hazard class: | (German Regulations) (Self-assessment) slightly hazardous for water. Do not allow undiluted product or large quantities of it to reach ground water, water course or sewer system. |
|---------------------|--|

**13. DISPOSAL CONSIDERATIONS:**

|  |   |
|--|---|
|  | <p>Recommendation: dispose of waste in accordance to applicable regional or local regulations.</p> <p>Uncleaned Packagings:<br/>Recommendation: Disposal must be made according to official regulations.</p> <p>Recommended cleansing agent: Water, if necessary with cleansing agents.</p> |
|--|---|

**14. Transport Information**

|                                       |                |  |  |
|---------------------------------------|----------------|--|--|
| DOT Regulations                       | Hazard Class   |  |  |
| Land transport ADR/RID (cross border) | ADR/RID class. |  |  |
| Maritime transport IMDG:              | IMDG class.    |  |  |
| Marine Pollutant:                     | No             |  |  |



|                                 |                                |
|---------------------------------|--------------------------------|
| <b>Commercial product name:</b> | <b>SED-CHEK®2 ESR CONTROLS</b> |
|---------------------------------|--------------------------------|

|                                   |  |
|-----------------------------------|--|
| <b>15. Regulatory information</b> |  |
| Hazard symbol:                    | T, I   |
| Warnings                          | See section 3: Hazards Identification  |
| Precautions:                      | Wear protective clothing upon contact, wash thoroughly with soap and water. Re-wash w/a detergent that contains anti microbial agents. |

|                              |   |
|------------------------------|---|
| <b>16. Other information</b> |   |
|                              | The information herein is believed to be correct as of the date hereof but is provided without warranty of any kind. The recipient of our products is responsible for observing any laws and guidelines applicable. |