

TECHNICAL DATA SHEET
Latex, Non-Sterile Medical Examination Glove
which meet ISO 11193-1 - Type 1

Intended Use

Non-sterile, single use, medical examination gloves are intended to prevent transmission of a wide variety of diseases to both patients and health care personnel. Specifically, gloves are intended to act as a barrier between the health care provider and patients to prevent contamination between patient and caregiver.

ISO 11193-1 - Type 1				
Characteristic			Acceptable Limit	
Dimension & Tolerance	Width (mm)	Extra-Small (XS)	≤ 80	
		Small (S)	70 - 90	
		Medium (M)	85 - 105	
		Large (L)	100 - 120	
		X-Large (XL)	≥ 110	
	Length (mm)	XS & S: ≥ 220 Others: ≥ 230		
Thickness (mm)	Palm	Smooth: 0,08-2,00	Textured: 0,11-2,03	
	Finger	Smooth: ≥ 0,08	Textured: ≥ 0,11	
Physical Requirement	Before ageing	Force at Break, N	≥ 7,0	
		Elongation at Break,%	≥ 650%	
	After ageing	Force at Break, N	≥ 6,0	
		Elongation at Break,%	≥ 500%	
Description			Ambidextrous. Refer to Appendix 1 for color and packaging information.	

Additional Information

Storage Condition	Country of Origin
Store in a cool, dry place. Avoid excessive heat (40°C, 104°F). Shield open box from direct sun and fluorescent light.	Malaysia

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1124	Medicom SafeTouch "Connect" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Chlorinated, Textured	Non-Sterile
1125	Medicom SafeTouch "Connect" Powdered Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Natural	Smooth	Non-Sterile
1162	Medicom SafeTouch "Connect" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Pink	Textured	Non-Sterile
1222	Medicom Safe Touch "UltraGrip" Powder-Free Latex Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case	A = X-small, B = Small, C = Medium, D = Large	Natural	Textured	Non-Sterile