

MODEL DHH20-30SDFMS

FLAMMABLE MATERIAL STORAGE FREEZER -18°C TO -25°C

Capacity [Cu. Ft.] 30 Cu. Ft.
Capacity [Liters] 849 Liters
Temperature Range [°C] -18°C to -25°C
Temperature Range [°F] 0°F to -13°F

Dimensions Width x Depth x Height
 Exterior [Inches] 35.75 x 35.5 x 72.5"
 Interior [Inches] 32" x 22.5" x 65"
 Exterior [Centimeters] 90 x 90 x 184
 Interior [Centimeters] 81 x 57 x 165

Compressor Hermetic Compressor
Air Circulation Gravity Flow
Defrost Method Manual
Refrigerant R404A
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 1 Door, Keyed Lock
Door Material Insulated Steel
Porthole 3/4" Diameter
Leveling Legs 4
Storage 4 Fixed Shelves
Shipping Weight 380 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 6A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
 Energy Star Rated

Controller

Display Type Mechanical Thermostat

» Digital Control Option

*Digital control option includes control & alarm system.
 Note: If a temperature recorder is purchased, the digital control option must also be purchased.

» Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor.