

MODEL C85-17X

EXPLOSION PROOF CHEST FREEZER -40°C TO -85°C

Capacity [Cu. Ft.]	17 Cu. Ft.
Capacity [Liters]	480 Liters
Temperature Range [°C]	-40°C to -85°C
Temperature Range [°F]	-40°F to -120°F
Dimensions	Width x Depth x Height
Exterior [Inches]	81.25" x 37.5" x 48.5"
Interior [Inches]	69" x 20" x 21"
Exterior [Centimeters]	206 x 95 x 123
Interior [Centimeters]	175 x 50 x 53
Refrigeration System	Cascade-Style System
Compressors	2, Hermetic Compressors
Compressor Container	Explosion Proof
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	HCFC & CFC Free
Exterior	Powder Coated Steel
Interior	Zinc Coated Galvanized Steel
Insulation	Polyurethane (Ecomate) Foam
Additional Insulation	Beaded Foam Sub-Lids
Lockable Lid	Lockable with Padlock
Casters	8, Swivel with Locks
Storage	Optional Inventory Racks
Shipping Weight	975 Lbs.

Electrical

Standard Voltage	115 Volts, 60 Hertz, 1 Phase
Amps / Run Amps / Plug (Optional Voltages)	20A 15A NEMA 5-20 208 Volts, 230 Volts
Power Configuration	Hardwired
Certification(s)	UL Rated
Certification File Number	SA10711
Certification Classification	Class I, Division I & II, Groups C & D

Controller

Display Type	Explosion Proof, Digital
--------------	--------------------------

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor