

MODEL DH4-12GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 9.19 Cu. Ft.
Capacity [Liters] 260 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions
Exterior [Inches] 25" x 22.75" x 63.25"
Interior [Inches] 22" x 16.5" x 43.75"
Exterior [Centimeters] 63 x 57 x 160
Interior [Centimeters] 55 x 41 x 111

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 1 Door, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter (Top Back Wall)
Casters 4
Storage 4 Adjustable Shelves
Shipping Weight 400 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 4.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed
Energy Star Rated

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.