

MODEL DH4-8GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	6.17 Cu. Ft.
Capacity [Liters]	175 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	19" x 20.75" x 71.5"
Interior [Inches]	15.75" x 14.5" x 46.75"
Exterior [Centimeters]	48 x 52 x 181
Interior [Centimeters]	40 x 36 x 118

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter (Top Back Wall)
Casters	4
Storage	5 Adjustable Shelves
Shipping Weight	300 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 2.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

QUOTED BY: