

MODEL DHK4-49GD-T

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 46.62 Cu. Ft.

Capacity [Liters] 1320 Liters

Temperature Range [°C] 2°C to 8°C

Temperature Range [°F] 35°F to 46°F

Dimensions Width x Depth x Height

Exterior [Inches] 57.25" x 34.75" x 78.75"

Interior [Inches] 52.75" x 26.25" x 59.5"

Exterior [Centimeters] 145 x 88 x 200

Interior [Centimeters] 134 x 67 x 152

Compressor Hermetic Compressor

Air Circulation Forced Air

Defrost Method Automatic

Refrigerant R290

Exterior Painted Steel

Interior Stainless Steel

Insulation US EPA and SNAP approved

Lockable Door Yes, Keyed

Porthole Yes, 1" Diameter

Casters 4, with Leveling Feet

Storage 10 Adjustable Shelves

Shipping Weight 800 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage Line | Running 15 Amp Dedicated | 10.06

Power Consumption 11.17 kW/d

NEMA Configuration NEMA 5-15

Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital

Touchscreen Yes, Color Display

Alarms High/Low Temp, Door Ajar

Power Failure, Low Battery

Sensor Failure, USB Failure

Data Logging Yes, Adjustable Intervals


Data Download Yes, via USB, PDF Format

Min / Max Temperature Yes, Display and Reset

Alarm Relay Yes, Dry Contacts

Backup Battery Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 2 Year Parts & Labor