

MODEL DHP4-38SGD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 38 Cu. Ft.
Capacity [Liters] 1076 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions
Exterior [Inches] 43.5" x 32" x 81"
Interior [Inches] 40.5" x 27" x 55"
Exterior [Centimeters] 110 x 81 x 205
Interior [Centimeters] 102 x 68 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 2 Sliding Doors, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter
Casters 4, Locking
Storage 10 Adjustable Shelves
Shipping Weight 550 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 5.0A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital
Touchscreen Yes, 4.3" Color Display
Alarms High/Low Temp
Alarm Notification Yes, Audible, Email, Text
Data Logging Yes, Adjustable Intervals
Data Download Yes, via USB, CSV Format
LAN Remote Access Yes, via VNC
Min / Max Temperature Yes, Display and Reset
Alarm Relay Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.