

## MODEL MV4-2UCRGDDA

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

**Capacity [ Cu. Ft. ]** 2.64 Cu. Ft.  
**Capacity [ Liters ]** 75 Liters  
**Temperature Range [ °C ]** 2°C to 8°C  
**Temperature Range [ °F ]** 35°F to 46°F

**Dimensions** Width x Depth x Height  
**Exterior [ Inches ]** 21.25" x 22" x 30"  
**Interior [ Inches ]** 17.5" x 17.5" x 21.75"  
**Exterior [ Centimeters ]** 54 x 56 x 76  
**Interior [ Centimeters ]** 44 x 44 x 54

**Compressor** Hermetic Compressor  
**Air Circulation** Forced Air  
**Defrost Method** Automatic  
**Refrigerant** R600a  
**Exterior** Steel  
**Interior** Polystyrene ( HIPS )  
**Insulation** US EPA and SNAP approved  
**Lockable Door** Yes, Keyed  
**Porthole** Yes, 1" Diameter  
**Casters** 2, Rear  
**Leveling Legs** 2, Front  
**Storage** 2 Adjustable Shelves, 1 Basket  
**Shipping Weight** 150 lbs.


#### Electrical

**Voltage:** 115 Volts 60 Hertz 1 Phase  
**Amperage Line | Running** 15 Amp Dedicated | 1.24A  
**Power Consumption** 1.24 kW/d  
**NEMA Configuration** NEMA 5-15  
**Certification(s)** UL Listed

#### Controller

**Display Type** Microprocessor, Digital  
**Alarms** High/Low Temp, Door Ajar  
Power Failure, Low Battery  
Sensor Failure, USB Failure  
**Data Logging** Yes, Adjustable Intervals  
**Data Download** Yes, via USB, PDF Format  
**Min / Max Temperature** Yes, Display and Reset  
**Alarm Relay** Yes, Dry Contacts  
**Backup Battery** Yes, Only Powers Alarm

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor