

Sheldon Forced Air Oven Model SMO28-2

**27.5 Cubic Ft. Upright
Temperature Range: To 260°C**

Specifications:

* Capacity	27.5 cu. Ft.
* Temperature Range	15°C above ambient to 260°C
* Dimensions	
* Exterior w x d x h	37.5 x 35.0 x 78.3" 95.3 x 88.9 x 198.8 (cm)
* Interior w x d x h	30.8 x 25.0 x 61.0" 78.2 x 63.5 x 157.5 (cm)

Design Information:

The SMO28-2 is the largest general purpose oven on the market (also consider the SMO38HP-2 high performance ovens for the best uniformity and functionality). The independent overtemperature safety feature is user adjustable and provides added safety. The turbo blower-assisted forced airflow design coupled with a heavy-duty blower motor offers excellent temperature uniformity and fast recovery times. Airflow is directed evenly across all shelves and samples to ensure consistent drying, curing and baking results.



Features:

Stainless Steel Shelves	Yes
3" Adjustable Exhaust Port	Yes
No Exposed Heating	Yes
Wrap Around Insulation	Yes
Digital Timer	Yes
Temperature Uniformity	3.0°C
Temperature Stability	0.2°C
Shelving	6 Supplied (16 Maximum)
Voltage	230/60/1
Amperage	20
Approval	CE Approved
Shipping Weight	495 lbs.
Warranty	Two Year

Specifications subject to change without notice