

Spirit®

HYDROCOLLOID ADHESIVE SHEATH

- First and only self-adhering male sheath to combine the excellent moisture-wicking properties of a hydrocolloid adhesive with the exceptional breathability of silicone^{1,2}
- A broad line of male external catheters with a diverse range of coverage
- Hydrocolloid adhesive wicks away moisture to help protect skin integrity, minimize skin maceration and maintain adhesion^{1,2}
- Made with 100% silicone



Style 1



Offers built-in adhesive in a length and placement found with traditional male external catheters.

Size		30/Box	100/Box
Small	25mm	35301	35101
Medium	29mm	35302	35102
Intermediate	32mm	35303	35103
Large	36mm	35304	35104
X-Large	41mm	35305	35105

Style 2



The 1.5" sheath with forward placement of adhesive designed to minimize unintended detachment.

Size		30/Box	100/Box
Small	25mm	37301	37101
Medium	29mm	37302	37102
Intermediate	32mm	37303	37103
Large	36mm	37304	37104
X-Large	41mm	37305	37105

Style 3



Our most popular style provides additional adhesive area compared to similar silicone self-adhering catheters. Combined with forward placement of adhesive designed to minimize unintended detachment.²

Size		30/Box	100/Box
Small	25mm	39301	39101
Medium	29mm	39302	39102
Intermediate	32mm	39303	39103
Large	36mm	39304	39104
X-Large	41mm	39305	39105