



MYCO
MEDICAL

TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET



TECHNOCUT® Disposable Scalpels

Technocut® Disposable Scalpels undergo 100% automated inspection to ensure consistent, high-quality scalpels with secure blade-to-handle attachment.

- Non-slip, clear shield improves safety and enables easy blade size identification
- Ergonomic, non-slip plastic handle with integrated metric ruler for convenience

- Scalpel blades produced from high-quality Swiss stainless steel
- 100% automated inspection verifies blade to handle attachment
- Cold-staking used for attaching blade to handle, improving consistency
- Scalpels manufactured in facilities that conform to ISO 13485

BRAND (Company)	Cold-Staking Blade-to-Handle Attachment	Non-Slip, Clear Shield	Ergonomic Handle	Full Range of Sizes
Technocut® Disposable Scalpels (MYCO Medical)	✓	✓	✓	✓
Bard-Parker™ Disposable Scalpels (Aspen)		✓	✓	
Disposable Scalpels (Miltex)		✓		



Ergonomic handle with integrated metric ruler



Sizes 10-25 available sterile



Cold staking used for blade-to-handle attachment

Sterile Technocut® Disposable Scalpels

BLADE SIZE	MYCO CODE	UNIT OF MEASURE
#10	6008T-10	10/bx, 100 bx/cs
#11	6008T-11	10/bx, 100 bx/cs
#12	6008T-12	10/bx, 100 bx/cs
#12B	6008T-12B	10/bx, 100 bx/cs
#15	6008T-15	10/bx, 100 bx/cs
#15C	6008T-15C*	10/bx, 100 bx/cs
#20	6008T-20	10/bx, 100 bx/cs
#21	6008T-21	10/bx, 100 bx/cs
#22	6008T-22	10/bx, 100 bx/cs
#23	6008T-23	10/bx, 100 bx/cs
#24	6008T-24*	10/bx, 100 bx/cs
#25	6008T-25*	10/bx, 100 bx/cs

*Special Order Item — Minimum orders and lead times may apply

Non-Sterile Technocut® Disposable Scalpels

BLADE SIZE	MYCO CODE	UNIT OF MEASURE
#10	6000CST-10	1000/case
#11	6000CST-11	1000/case
#15	6000CST-15	1000/case
#20	6000CST-20	1000/case
#21	6000CST-21	1000/case
#22	6000CST-22	1000/case

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



MYCO
MEDICAL
TheSharpChoice®

2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com

www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017, 2024 MYCO Medical, Inc. All rights reserved.
MYCO® is a registered trademark of MYCO Medical, Inc.

LIT-00228 Rev 2 Eff 06.19.2024