

Thumper® Equine Pro



Based on the design and features of the Maxi Pro, the Equine Professional Massager delivers even more power and strength, so it can be used on horses. The Equine Pro is the only massager capable of improving your horse's overall performance.

Thumper has been used with remarkable success on race horses, hunter/jumpers, rodeo and reining horses, dressage horses, eventing, professional cutting horses, polo ponies, draft horses and chuck wagon teams.

Thumper relieves muscle aches and improves the overall attitude of physically strained horses. The horse will enjoy being massaged, and both horse and owner will benefit from the results.

Features

- 12 ft. (3.66 m) Power Cord
- Variable Speed: 18 - 35 Pulses/ Sec
- Stroke: 0.32" (8.1 mm)
- Class 1 FDA Medical Device
- Ground Fault Current Interrupter Device Included

Available Accessories

- Carry Case
- Foot Cushion
- Maxi Pro Wall Mount



	Weight		Dimensions	
Product	7 lbs	3.2 kg	10 x 5.5 x 7.5"	25.4 x 14.0 x 19.1 cm
Single Pack	9 lbs	4.1 kg	12 x 9 x 10"	30.5 x 22.9 x 25.4 cm
Case Pack (4 units)	36 lbs	16.3 kg	12 x 8 x 10"	30.5 x 45.7 x 50.8 cm

	Model Number	UPC Number	Electrical
North America (NA)	1007ENA	779665002512	120 Volts, 60 Hertz
Europe (EU)	1007EEU	779665002529	230 Volts, 50 Hertz
United Kingdom (UK)	1007EUK	779665002536	230 Volts, 50 Hertz
Australia (AU)	1007EAU	779665002543	240 Volts, 50 Hertz

forward all orders to : orderdesk@thumpermassager.com



1-800-THUMPER (848-6737)
1-905-477-5222 (Int'l)
1-905-477-5329 (Fax)
info@thumpermassager.com
www.thumpermassager.com