

Thumper® Lithium8

The Lithium8 is built for professional use and it's ideal for chiropractors, physiotherapists and sports massage therapists. Great to use at your office, clinic, football pitch, or the ski chalet. The Lithium8 can provide a full body massage in 5 to 8 minutes, with no need of a wall outlet or no chance of a power cable getting in the way of your treatment.

With its unique dual-drive assembly system and the new lithium battery technology Lithium8 is just as strong and effective as plug-in models. The balanced design and the weight of 7.9 pounds eliminate the need to apply pressure and reduces feedback to the hand for effortless handling. With a depth of stroke of 6mm it delivers a very strong massage.

The batteries last approximately 60 minutes, about 10 massage sessions. The easy swappable battery allows for continuous use. The battery has a quick release and charges on a separate charging nest.



Features

- Package includes Massager, Charging Nest, and 1 Lithium Battery
- 5 Speeds: 20 - 35 Pulses/ Sec
- Stroke: 0.24" (6.0 mm)
- Class 1 FDA Medical Device

Available Accessories

- Thumper Lithium Battery – so you can use one while the other charges
- Lithium8 Massager Carry Case – to take it with you
- Additional Charging Nest

	Weight		Dimensions	
Product	7.9 lbs	3.4 kg	10.6x10.2x7.9"	27x26x20 cm
Single Pack	10.4 lbs	4.7 kg	13.8x11x8.3"	35x28x21 cm
Case Pack (4 units)	42 lbs	18.8 kg	21.7x13.9x17"	55x35x43 cm

	Model Number	UPC Number
Lithium8 package (Batt + Charger)	J501	7 79665 00801 9
Battery only	J300	7 79665 00803 3
Charging Nest only	J60	7 79665 00802 6

forward all orders to : orderdesk@thumpermassager.com



1-800-THUMPER (848-6737)
1-905-477-5222 (Int'l)
1-905-477-5329 (Fax)
info@thumpermassager.com
www.thumpermassager.com