

# Thumper® Maxi Pro



This is the authentic deep tissue massager, originally designed for clinical use. The Maxi Pro is ideal for warming up muscles prior to Chiropractic adjustment. The Maxi Pro's 8 massage spheres cover a broad area, ideal for large muscle groups. A full-body massage takes only 5-8 minutes.

Maxi Pro is designed to direct 95% of its energy into the patient's muscles, assuring virtually no feedback on the therapist's hands. No need for oils or lotions. Thumper® works through clothing, providing flexibility for use in all clinical and professional environments.

This product offers limited use for self-massage.

## Features

- 12 ft. (3.66 m) Power Cord
- Variable Speed: 18 - 35 Pulses/ Sec
- Stroke: 0.25" (6.4 mm)
- Class 1 FDA Medical Device

## Available Accessories

- Carry Case
- Foot Cushion
- Maxi Pro Wall Mount



	Weight		Dimensions	
<b>Product</b>	7 lbs	3.2 kg	10 x 5.5 x 7.5"	25.4 x 13.9 x 19.1 cm
<b>Single Pack</b>	9 lbs	4.1 kg	12 x 9 x 10"	30.5 x 22.9 x 25.4 cm
<b>Case Pack (4 units)</b>	36 lbs	16.3 kg	12 x 18 x 20"	30.5 x 45.7 x 50.8 cm

	Model Number	UPC Number	Electrical
<b>North America (NA)</b>	1007NA	779665002000	120 Volts, 60 Hertz
<b>Europe (EU)</b>	1007EU	779665002024	230 Volts, 50 Hertz
<b>United Kingdom (UK)</b>	1007UK	779665002031	230 Volts, 50 Hertz
<b>Australia (AU)</b>	1007AU	779665002048	240 Volts, 50 Hertz

forward all orders to : [orderdesk@thumpermassager.com](mailto:orderdesk@thumpermassager.com)



**1-800-THUMPER (848-6737)**  
**1-905-477-5222 (Int'l)**  
**1-905-477-5329 (Fax)**  
[info@thumpermassager.com](mailto:info@thumpermassager.com)  
[www.thumpermassager.com](http://www.thumpermassager.com)