

# Thumper® Mini Pro

The Mini Pro is a Professional Grade product designed for durability and effectiveness. It is used around the world in both clinical and home environments. A strong massage with no work: Thumper Mini Pro does the work, not the user.

Its long handle design is perfect for treating others and for self-use at home. You can reach areas like the shoulders, hamstrings and mid-back. This product is excellent for self-care.

With three electronically controlled speed settings at 20, 30 and 40 pulses per second, the Mini Pro can be adjusted to accommodate personal preference. The Mini Pro is easy to hold and maneuver and over 95% of the units output is directed through the massage spheres eliminating feedback through the handle.



## Features

- 12 ft. (3.66 m) Power Cord
- Variable Speed: 20 - 40 Pulses/ Sec
- Stroke: 0.24" (6.0 mm)
- Class 1 FDA Medical Device

## Available Accessories

- Carry Case
- Mini Pro Wall Mount

	Weight		Dimensions	
<b>Product</b>	3 lbs	1.4 kg	13.5 x 6.5 x 5.25"	34.3 x 15.9 x 13.3 cm
<b>Single Pack</b>	5 lbs	2.3 kg	17 x 7 x 6"	43.2 x 17.8 x 15.4 cm
<b>Case Pack (6 units)</b>	30 lbs	13.6 kg	17 x 14 x 17"	43.2 x 35.6 x 43.2 cm

	Model Number	UPC Number	Electrical
<b>North America (NA)</b>	D401NA	779665001218	120 Volts, 60 Hertz
<b>Europe (EU)</b>	D401EU	779665001225	230 Volts, 50 Hertz
<b>United Kingdom (UK)</b>	D401UK	779665001232	230 Volts, 50 Hertz
<b>Australia (AU)</b>	D401AU	779665001249	240 Volts, 50 Hertz

forward all orders to : [orderdesk@thumpermassager.com](mailto:orderdesk@thumpermassager.com)



**1-800-THUMPER (848-6737)**  
**1-905-477-5222 (Int'l)**  
**1-905-477-5329 (Fax)**  
[info@thumpermassager.com](mailto:info@thumpermassager.com)  
[www.thumpermassager.com](http://www.thumpermassager.com)