

Thumper® Sport



The Sport is a cost effective, high quality massager for home-use. Weighing only 3 pounds it's easy to carry, and can be used at home, at the gym, or anywhere you travel. Its long ergonomic handle allows you to massage all areas of the body with ease.

The proven Thumper® drive system delivers a comfortable and penetrating tapotement action, with virtually no energy kick back to the handle. This makes the massager easy to hold, maneuver and operate.

Customize the massage experience through the use of soft or hard massage spheres. Use the variable speed control to select the most comfortable speed setting between 20 and 40 pulses per second.

Features

- 12 ft. (3.66 m) Power Cord
- Variable Speed: 20 - 40 Pulses/ Sec
- Stroke: 0.2" (5 mm)
- Class 1 FDA Medical Device
- 1 Set of Interchangeable Massage Skins

Available Accessories

- Carry Case
- Thumper Sport Wall Mount



	Weight		Dimensions	
Product	3 lbs	1.4 kg	6 x 4 x 15"	15.2 x 10.2 x 38.1 cm
Single Product	5 lbs	2.3 kg	7 x 6 x 20"	17.8 x 15.2 x 50.8 cm
Case Pack (6 units)	30 lbs	13.6 kg	20 x 11 x 20"	50.8 x 27.9 x 50.8 cm

	Model Number	UPC Number	Electrical
North America (NA)	E501NA	779665005018	120 Volts, 60 Hertz
Europe (EU)	E501EU	779665005025	230 Volts, 50 Hertz
United Kingdom (UK)	E501UK	779665005032	230 Volts, 50 Hertz
Australia (AU)	E501AU	779665005049	240 Volts, 50 Hertz

forward all orders to : orderdesk@thumpermassager.com



1-800-THUMPER (848-6737)
1-905-477-5222 (Int'l)
1-905-477-5329 (Fax)
info@thumpermassager.com
www.thumpermassager.com