

Thumper® Verve

A perfect choice for the home user, now available at cost-effective pricing yet built to professional quality standards.

Thumper® Verve is a single sphere massager that uses Thumper's proven percussive action. It is the smallest massager in Thumper's product line, fitting right into the palm of your hand and weighing only 2 lbs.

Thumper® Verve is ideal for treating precise, hard-to-get-at areas, and trigger points. You can target your massage on tight muscles and knots, releasing built-up tension in seconds.

Thumper® Verve is easy to use. Think of it as an extension of your hand, as if you were wearing a mitten. Simply put your hand over it, as if you were using a computer mouse, and strap the flexible band across the back of your hand. You don't need to grip it; it just stays on your hand. Glide the massager over your skin or clothes on its soft and smooth heel. Guide the massage head to the area you want to massage and let it do its work. For a deeper massage, simply apply some pressure with your fingertips. When used properly, there is virtually no feedback onto your hand or wrist.

Thumper® Verve is very ergonomic and intuitive; with a bit of practice you will be able to release the tightness from your muscles and feel relaxed and invigorated!

Features

- 12 ft. (3.66 m) Power Cord
- Variable Speed: 20 - 40 Pulses/ Sec
- Stroke: 0.20" (5.0 mm)
- Class 1 FDA Medical Device
- 3 Interchangeable Massage Skins
- Carry Bag Included

Specifications

- Size: 3.5"W x 7.5"L x 3.5"H (9cm x 19cm x 9cm).
- Weight: 2 lbs. (0.9 kg)
- Electrically approved UL, CSA, CE
- 2-year warranty



	Weight		Dimensions	
Product	2 lbs	0.9 kg	3.5 x 7.5 x 3.5"	8.9 x 19.1 x 8.9 cm
Single Pack	3 lbs	1.4 kg	7 x 6 x 11"	17.8 x 15.2 x 27.9 cm
Case Pack (6 units)	20 lbs	9.1 kg	20 x 11 x 10"	50.8 x 27.9 x 25.4 cm

	Model Number	UPC Number	Electrical
North America (NA)	G501NA	779665007012	120 Volts, 60 Hertz
Europe (EU)	G501EU	779665007029	230 Volts, 50 Hertz
United Kingdom (UK)	G501UK	779665007036	230 Volts, 50 Hertz
Australia (AU)	G501AU	779665007043	240 Volts, 50 Hertz

forward all orders to : orderdesk@thumpermassager.com



1-800-THUMPER (848-6737)
1-905-477-5222 (Int'l)
1-905-477-5329 (Fax)
info@thumpermassager.com
www.thumpermassager.com