



Specifications

Sky Blue Midnight Blue	Regular 21" (54 cm) - 600-21000X	Low 19" (48 cm) - 600-19000X
Overall width	24.75" (63 cm)	
Folded width	10.5" (26.7 cm)	
Length	28.25" (71.75 cm)	
Handle height range	32.4" - 40.2" (82.3 cm - 102 cm)	29.4" - 37.2" (74.7 cm - 95 cm)
Distance between handles	17.75" (45 cm)	
Seat height	21" (54 cm)	19" (48 cm)
Seat width	16" (40.6 cm)	
Seat depth (rollator)	12.5" (30.5 cm)	
Seat depth (transport chair)	16" (40.6 cm)	
Weight (rollator)	30 lbs (13.6 kg)	
Weight (total)	32 lbs (14.5 kg)	
Weight Capacity	275 lbs (125 kg)	