

# 2025

Viscot is proud to be a veteran owned small business (VOSB)

# VISCOT<sup>®</sup>

MEDICAL, LLC.

**The *Leader* In Surgical Marker Innovation**

Surgical Markers • Medication Labels • Neutral Zones • Drapes • Urinals

[www.viscot.com](http://www.viscot.com)

P: 1.800.221.0658 | P: 973.221.0658 | F: 973.887.3961 | Email: [viscotcs@viscot.com](mailto:viscotcs@viscot.com)

## ABOUT US

Viscot Medical, LLC, is the global leader in surgical site marking innovations. Viscot is committed to providing the healthcare industry with high quality disposable products at an economical price. We are proud to serve our global customers as a U.S. manufacturer for over 50 years.

Viscot Medical is dedicated to listening to clinicians and coming up with site marking solutions to meet each clinical need.

Viscot Medical is a veteran-owned, small business.

ISO 13485 CERTIFIED  
MDR CERTIFIED



## TABLE OF CONTENTS


About Our Inks .....	4
About Our Tips .....	5
About Our Barrels .....	6
New Products .....	8
Traditional Ink Surgical Markers .....	9
XL Prep Resistant Ink Markers .....	15
Non-Surgical Aesthetic Skin Markers .....	19
Specialized Site Markers .....	21
Utility Markers .....	24
Pre-Printed Medication Labels .....	27
Sharps & Staff Safety .....	34
Additional Products .....	36
How To Read Your Packaging .....	41

# SAME MARKERS. DIFFERENT PACKAGE.

Viscot Medical, LLC, is committed to creating a sustainable global future.



Over 20 years ago, Viscot introduced mini markers with a core goal of reducing waste. To support our goal, we are replacing our plastic tub with an environmentally friendly biodegradable box. Same product quality in a new eco-friendly package.



**TRADITIONAL INK**

Traditional ink is Gentian Violet ink and has been used in the industry to mark patients for over 40 years. It can withstand betadine and iodine scrubs. AllSkin<sup>®</sup> (pg 9) features a brilliant blue option available for visibility on all skin tones.

*Ideal for procedures above the neck or where betadine and iodine are used.*

*Featured on pages 9-13*



**GVFREE<sup>®</sup> INK**

Viscot Medical's GVFree<sup>®</sup> ink is a gentian violet free alternative to Viscot's Traditional ink. Eliminates concerns over Gentian Violet while maintaining the classic purple color.

*Ideal for procedures above the neck or where betadine and iodine are used.*

*Featured on page 14*



**ViXL<sup>®</sup> PREP RESISTANT INK**

Viscot's Patented ViXL<sup>®</sup> Prep Resistant Ink is visible up to 10x longer vs. traditional gentian violet ink. Able to withstand a 60 second prep of Chloraprep, ViXL<sup>®</sup> Ink remains visible against all common alcohol based prep solutions. Available in gentian violet free red and green ink for enhanced visibility on all skin tones.

*Ideal for any procedure where aggressive preps are used.*

*Featured on pages 15-18*



**EZ REMOVABLE INK<sup>®</sup>**

Viscot Medical's EZ Removable Ink is designed to come off easily with an alcohol swab or even warm soap and water. Available in red, green, and white.

*Ideal for non-surgical aesthetic procedures.*

*Featured on page 20*



**WATERPROOF UTILITY INK**

Black Utility ink is designed for writing on medication labels. Smudge-proof, waterproof, and fluid resistant.

*Ideal for writing on labels, glass, vinyl and other surfaces.*

*Featured on page 23-25*



**ULTRA FINE TIPS**

Viscot Medical's soft rounded Ultra Fine tips create fine lines while eliminating concerns over causing trauma to delicate tissues and reducing frustration from sharp marker points snagging making it difficult to draw.

*Ideal for facial procedures such as blepharoplasty procedures.*



**REGULAR TIPS**

This is the most common tip style found on the market. Designed to provide versatility and a middle-ground between fine lines and thick lines.

*Ideal for most general procedures and for correct site marking.*



**BOLD TIPS**

Viscot Medical's bold tips cover a wide surface area and create broad lines. Additional tipsurface makes it more difficult for the marker to clog and appear dry.

*Ideal for abdominoplasty procedures and vascular procedures.*



**BROAD TIPS**

Viscot Medical's broad tip is unique to our 1404 and 1405 markers. This tip was specifically designed to replace gentian violet ink stamp pads. The paintbrush like tip coats instruments easily and completely.

*Ideal for coating surgical Instruments with ink but especially cataract procedures.*

## ABOUT OUR BARRELS



### MINI MARKERS

Viscot Mini Markers are economical and eco-friendly. They're about half the size of a standard marker so they're about half the cost and produce half the waste.

*Perfect for promoting single patient use.*



### VALUE BARREL

This barrel feels like a pencil when gripped. Value barrels are thinner and shorter than Precision and Wide Barrels, but not as small as a Mini Marker. Less plastic means more cost savings.



### MULTI-TIP BARRELS

Features two tips one on each end. Multi-Tip barrels differ depending on the exact product numbers but are an economical, convenient and eco-friendly alternative by providing two markers in one.

Usually featuring different tips, the unique Blephmarker<sup>®</sup> (page 11) features two ultra fine tips, offering a fresh tip per eyelid.



### PRECISION BARREL

Our standard barrels designed with a comfortable grip in mind. Precision Barrels feel and perform like a fountain pen. Featured on our Precision markers and full size XL markers.



### WIDE BARREL

Thicker barrels that feel more like a standard marker and less like a pen. All Wide Barrel Markers have a clip cap.



### PRESSURIZED BARREL

Found in Viscot Medical's Visflow Markers. Designed to keep marking even through fluids on wet tissue.

*Ideal for abdominoplasty procedures, breast procedures, and liposuction.*

13.4 x 1 cm



**Precision/Classic Barrel**

- 1400 / 1401 - GV Ink | Reg Tip | NS, SR, SRL
- 1402 / 1403 - GV Ink | UF Tip | NS, SR
- 1404 / 1405 - GV Ink | Broad Tip | NS, SR
- 1406 / 1407 - Blue Ink | Reg Tip | NS, SR
- 1447XL - XL, XLR, XLG Ink | Reg Tip | NS, NSRL, SR, SRL
- 1448XL - XL Ink | UF Tip | NS, SR
- 1449XL - XL Ink | Bold Tip | NS, SR
- Visi-Tip - GV Ink | Radiopaque Tip | NS

13.6 x 1.1 cm



**Wide Barrel**

- 1408 - Perm Ink | Reg Tip | NS, SRL
- 1411 - Perm Ink | Fine Tip | NS, S, SRL
- 1416 - Perm Ink | Super Fine Tip | NS, S
- 1460 - GV Ink | Reg Tip | NS, SRL
- 1480 - GV Ink | Reg Tip | NS

12.9 x 0.5 cm



**Value Barrel**

- 1436 - GV Ink | UF Tip | NS, SR, SRL
- 1437 - GV, GVF, XLR Ink | Reg Tip | NS, NSR, NSRL, SR, SRL

13.4 x 0.8 cm



**1438 Barrel**

- GV Ink | Reg Tip | NS, SR, SRL



**Pressurized Barrel - Squeeze**

- VF401 - GV Ink | Bold Tip | SR
- VF500 - GV Ink | Reg Tip | SR
- VF600 - GV Ink | UF Tip | SR

10 x 0.5 cm



**Rounded Mini Barrel**

- 1450XL - XL, XLR Ink | Reg Tip | NS, NSW, SR
- 1456XL - XL Ink | UF Tip | NS, NSW, SR
- 1481GVF - GVF Ink | Reg Tip | NS, SR



**Hexagon Mini Barrel**

- 1451 - GV Ink | Reg Tip | NS, NSW, SR

10 x 0.5 cm



**Wide Mini Barrel**

- 1440 - EZW Ink | Reg Tip | NS
- 1444 - EZG Ink | Reg Tip | NS
- 1445 - EZR Ink | Reg Tip | NS
- 1453 - Perm Ink | Reg Tip | NS, SRL

14.2 x 0.5 cm



**1424 Barrel**

- GV Ink | 2 UF Tips | NS, SR

13.5 x 0.7 cm



**1425 Barrel**

- GV Ink | Reg & UF Tips | NS, NSRL, SR

13.5 x 0.7 cm



**1422 Barrel**

- GV & Perm Ink | Regular & Fine Tips | NS, NSR, NSRL, S, SR, SRL



**1462 Barrel**

- XL & Perm Ink | Regular & Fine Tips | SR, SRL



## NEW PRODUCTS

### WM200 WOUND MARKER WITH GENTIAN VIOLET FREE INK Regular Tip | Mini Barrel



PART #	DESCRIPTION	CASE QTY
WM200-200	Wound Marker: skin marker designed for marking and monitoring infection sites, wound margins, etc. Regular tip purple gentian violet ink marker.	200

### VF401 VISFLOW ULTRA SQUEEZABLE INK MARKER Bold Tip | Pressurized Barrel | Traditional Ink



PART #	DESCRIPTION	CASE QTY
VF401SR-50	Sterile Bold Tip w/ Ruler	50
VF401SR-100	Sterile Bold Fine Tip w/ Ruler - 2/Bx of 50 Markers	100



## **TRADITIONAL INK SURGICAL MARKERS**

Ideal for Pre-Operative procedures where Betadine or Iodine are used.  
All Sterile versions are ideal for Intraoperative use.



# TRADITIONAL INK SURGICAL MARKERS



## 1400/1401 PRECISION SURGICAL MARKER Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1400-100	Sterile Regular Tip w/ Ruler	100
1400-100L	Sterile Regular Tip w/ 4 Labels & Ruler	100
1401-100	Non-Sterile Regular Tip	100



## 1402/1403 PRECISION SURGICAL MARKER Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1402-100	Sterile Ultra Fine Tip w/ Ruler	100
1403-100	Non-Sterile Ultra Fine Tip	100



## 1404/1405 PRECISION SURGICAL MARKER Broad Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1404-100	Sterile Broad Tip w/ Ruler	100
1405-100	Non-Sterile Broad Tip	100



## 1406/1407 PRECISION SURGICAL MARKER / ALLSKIN Brilliant Blue Ink | Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1406-100	Non-Sterile Regular Tip	100
1407-100	Sterile Regular Tip w/ Ruler	100



## 1422 DUAL TIP BARREL Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1422-100	Non-Sterile 1 Regular & 1 Fine Tip	100
1422S-100	Sterile 1 Regular & 1 Fine Tip	100
1422SR-100	Sterile w/ Ruler 1 Regular & 1 Fine Tip	100
1422SRL9-100	Sterile 1 Regular & 1 Fine Tip w/ Ruler & 9 Blank Labels	100



# TRADITIONAL INK SURGICAL MARKERS



## 1424 TWIN TIP BARREL

Blephmarker: 2 Ultra Fine Tips | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1424-100	Non-Sterile 2 Ultra Fine Tips	100
1424SR-100	Sterile 2 Ultra Fine Tips w/ Ruler	100



## 1425 TWIN TIP BARREL

Regular Tip & Ultra Fine Tip | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1425-100	Non Sterile 1 Regular & 1 Ultra Fine Tip	100
1425SR-100	Sterile 1 Regular & 1 Ultra Fine Tip w/ Ruler	100
1425SRL9-100	Sterile 1 Regular & 1 Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



## 1436 VALUE SURGICAL MARKER

Ultra Fine Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1436-100	Non Sterile Ultra Fine Tip	100
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	100
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



## 1437 VALUE SURGICAL MARKER

Regular Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1437-100	Non-Sterile Regular Tip	100
1437SR-100	Sterile Regular Tip w/ Ruler	100
1437SRL6-100	Sterile Regular Tip w/ Ruler & 6 Blank Labels	100
1437SRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



## 1438 VALUE SURGICAL MARKER

Regular Tip | Value Barrel | Clip Cap

PART #	DESCRIPTION	CASE QTY
1438SR-100	Sterile Regular Tip, Clip Cap w/ Ruler	100
1438SRL9-100	Sterile Regular Tip, Clip Cap w/ Ruler & 9 Labels	100



# TRADITIONAL INK SURGICAL MARKERS



## 1451 MINI SURGICAL MARKER Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1451-200	Non-Sterile Regular Tip	200
1451-1000	Non-Sterile Regular Tip	1,000
1451SR-100	Sterile Regular Tip w/ Ruler	100
1451W-100	Non-Sterile Regular Tip Individually Wrapped (Proof of Use)	100



## 1460 STANDARD TRADITIONAL INK SURGICAL MARKER Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1460-100	Non-Sterile Regular Tip	100
1460SRL-100	Sterile Regular Tip w/ Ruler & Labels & 4 Labels	100
1460SRL9-100	Sterile Regular Tip w/ Ruler & 9 Labels	100




## 1480 STANDARD SURGICAL MARKER Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1480-100	Non-Sterile Regular Tip	100



## MMF MINI SURGICAL MARKER Ultra Fine Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
MMF050C	Non-Sterile Ultra Fine Tip	300
MMF1000B	Non-Sterile Ultra Fine Tip	1,000



## MMX MINI MAX SURGICAL MARKER Regular Tip | Mini Max Barrel

PART #	DESCRIPTION	CASE QTY
MMX250	Non-Sterile Regular Tip Individually Wrapped (Proof of Use)	500

# TRADITIONAL INK SURGICAL MARKERS




## MMX MINI MAX SURGICAL MARKER (Continued)

<b>MMX1000</b>	Non-Sterile Regular Tip Individually Wrapped (Proof of Use)	1,000
<b>MMX100C</b>	Non-Sterile Regular Tip	400
<b>MMX1000B</b>	Non-Sterile Regular Tip	1,000
<b>MMX2000B</b>	Non-Sterile Regular Tip	2,000
<b>PSX200</b>	Sterile Regular Tip	200
<b>PSX200R</b>	Sterile Regular Tip w/ Ruler	200



## VF401 VISFLOW ULTRA SQUEEZABLE INK MARKER Bold Tip | Pressurized Barrel

<b>VF401SR-50</b>	Sterile Bold Tip w/ Ruler	50
<b>VF401SR-100</b>	Sterile Bold Fine Tip w/ Ruler - 2/Bx of 50 Markers	100



## VF500 VISFLOW ULTRA SQUEEZABLE INK MARKER Regular Tip | Pressurized Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>VF500SR-50</b>	Sterile Regular Tip w/ Ruler	50
<b>VF500SR-100</b>	Sterile Regular Tip w/ Ruler - 2/Bx of 50 Markers	100



## VF600 VISFLOW ULTRA SQUEEZABLE INK MARKER Ultra Fine Tip | Pressurized Barrel

<b>PART #</b>	<b>DESCRIPTION</b>	<b>CASE QTY</b>
<b>VF600SR-50</b>	Sterile Ultra Fine Tip w/ Ruler	50
<b>VF600SR-100</b>	Sterile Ultra Fine Tip w/ Ruler - 2/Bx of 50 Markers	100

# GENTIAN VIOLET FREE TRADITIONAL INK MARKERS

A Gentian Violet Free Alternative To Traditional Gentian Violet Free Markers. Complies with California Prop 65.



**1437 VALUE GVFREE® INK SKIN MARKER**  
Regular Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1437GVFSR-100	Sterile Regular Tip w/ Ruler	100
1437GVFSRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



**1481 MINI GVFREE® INK SKIN MARKER**  
Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1481GVF-200	Non-Sterile Regular Tip	200
1481GVF-1000	Non-Sterile Regular Tip	1,000
1481GVFSR-100	Sterile Regular Tip w/ Ruler	100



## **XL PREP RESISTANT INK MARKERS**

Patented Prep Resistant Ink. Designed For Pre-Operative Site Marking Where Aggressive Preps Are Used. Visible Throughout Prepping with All Prep Including Chlorhexidine & Alcohol Based Preps.\*



\* Products with XL ink or ViXL® covered by US Patents 11203697 11359103 and other patents pending.

# XL PREP RESISTANT INK MARKERS



## 1437 XL PREP RESISTANT INK VALUE SKIN MARKER

Regular Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1437XLSR-100	Sterile Regular Tip w/ Ruler	100
1437XLSRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



## 1447 XL PREP RESISTANT INK SKIN MARKER

Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1447XL-50	Non-Sterile Regular Tip	50
1447XLMC-50	Non-Sterile Regular Tip: 18 Red, 16 Green., & 16 Purple	50
1447XLSR-100	Sterile Regular Tip w/ Ruler	100



## 1448 XL PREP RESISTANT INK SKIN MARKER

Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1448XL-50	Non-Sterile Ultra Fine Tip	50
1448XLSR-100	Sterile Ultra Fine Tip w/ Ruler	100



## 1449 XL PRE-SURGICAL MARKER

Bold Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1449XL-50	Non-Sterile Bold Tip	50
1449XLSR-100	Sterile Bold Tip w/ Ruler	100



## 1450 MINI MARKER XL PREP RESISTANT INK SKIN MARKER

Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1450XL-200	Non-Sterile Regular Tip	200
1450XL-1000	Non-Sterile Regular Tip	1,000
1450XLSR-100	Sterile Regular Tip w/ Ruler	100
1450XLW-100	Non-Sterile Regular Tip Individually Wrapped (Proof of Use)	100



# XL PREP RESISTANT INK MARKERS



## 1456 MINI MARKER XL PREP RESISTANT INK SKIN MARKER

Ultra Fine Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1456XL-200	Non-Sterile Ultra Fine Tip	200
1456XL-1000	Non-Sterile Ultra Fine Tip	1,000
1456XLSR-100	Sterile Ultra Fine Tip w/ Ruler	100
1456XLW-100	Non-Sterile Ultra Fine Tip Individually Wrapped (Proof of Use)	100



## 1447 XL PREP RESISTANT RED INK SKIN MARKER

Regular Tip | Precision Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1447XLR-50	Non-Sterile Ultra Fine Tip	50
1447XLMC-50	Non-Sterile Regular Tip	50



## 1447 XL PREP RESISTANT GREEN INK SKIN MARKER

Regular Tip | Precision Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1447XLG-50	Non-Sterile Regular Tip	50
1447XLMC-50	Non-Sterile Regular Tip	50
1447XLGSR-100	Sterile Regular Tip w/ Ruler	100



## 1462 DUAL TIP BARREL XL PREP RESISTANT INK SKIN MARKER

Regular Tip + XL Prep Resistant Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1462XLSR-100	Sterile w/ Ruler 1 Regular & 1 Fine Tip	100
1462XLSRL9-100	Sterile 1 Regular & 1 Fine Tip w/ Ruler & 9 Blank Labels	100

# XL PREP RESISTANT INK MARKERS



## 1450 XL PREP RESISTANT RED INK MINI SKIN MARKER

Regular Tip | Mini Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1450XLRED-1000	Non-Sterile Regular Tip	1,000
1450XLRED-200	Non-Sterile Regular Tip	200
1450XLREDSR-100	Sterile Regular Tip w/ Ruler	100



## 1437 XL PREP RESISTANT RED INK VALUE SKIN MARKER

Regular Tip | Value Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1437XLREDSR-100	Sterile Regular Tip w/ Ruler	100
1437XLREDSRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100





## **NON-SURGICAL AESTHETIC SKIN MARKERS**

Marking Non-Surgical Aesthetic Procedures  
Like Consultations, Botox, & Fillers.



# NON-SURGICAL AESTHETIC SKIN MARKERS



## 1440 EZ REMOVABLE WHITE INK SKIN MINI MARKER

Regular Tip | Mini Barrel | White Ink

PART #	DESCRIPTION
1440-100	Non-Sterile Regular Tip
1440-500	Non-Sterile Regular Tip

CASE QTY



100  
500



## 1444 EZ REMOVABLE GREEN INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Green Ink

PART #	DESCRIPTION
1444-100	Non-Sterile Regular Tip
1444-500	Non-Sterile Regular Tip

CASE QTY



100  
500



## 1445 EZ REMOVABLE RED INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Red Ink

PART #	DESCRIPTION
1445-100	Non-Sterile Regular Tip
1445-500	Non-Sterile Regular Tip

CASE QTY



100  
500



# **SPECIALIZED SITE MARKERS**

Alternatives To Surgical Markers.



# SPECIALIZED SITE MARKERS



## AM106 ACCU-MARK RADIOPAQUE SITE MARKER

Marking Under Fluoroscopy

PART #	DESCRIPTION	CASE QTY
AM106-50	Non-Sterile w/ Traditional Surgical Ink	50
AM106S-50	Sterile w/ Traditional Surgical Ink	50
AM106-100	Non-Sterile w/ Traditional Surgical Ink	100



## SMALL 1/2" YES SITE MARKERS

Cataract, Pediatric & Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
TYS-100W	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	100
TYS-400	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Per Box / 4Box / 1 Cs	400



## LARGE 1" YES SITE MARKERS

Cataract, Pediatric & Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
TYL-100W	Non-Sterile Temporary 1" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	100
TYL-400	Non-Sterile Temporary 1" in Green "Yes" - 100 Per Box / 4 Box / 1 Cs	400

# SPECIALIZED SITE MARKERS

**NEW**

## **TOR-MARK-V**

Disposable Cataract Marking Kit

Mini Max Regular Tip Traditional Ink Marker, Dial with Level and Indexed Cap 5/cs





# UTILITY MARKERS



## 1408 WATERPROOF PERMANENT MARKER

Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1408-50	Non-Sterile Regular Tip	50
1408SRL-100	Sterile Regular Tip w/ Ruler & 4 Blank Labels	100



## 1411 WATERPROOF PERMANENT MARKER

Fine Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1411-100	Non-Sterile Fine Tip	100
1411S-100	Sterile Fine Tip	100
1411SRL9-100	Sterile Fine Tip w/ Ruler & 9 Blank Labels	100



## 1416 WATERPROOF PERMANENT MARKER

Ultra Fine Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1416-100	Non-Sterile Superfine Tip	100
1416S-100	Sterile Superfine Tip	100



## 1422 DUAL TIP SURGICAL INK MARKER & PERMANENT INK MARKER

Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1422-100	Non-Sterile Regular Tip & Non-Sterile Fine Tip	100
1422S-100	Sterile Regular Tip & Sterile Fine Tip	100
1422SR-100	Sterile Regular Tip & Sterile Fine Tip w/ Ruler	100
1422SRL9-100	Sterile Regular Tip & Sterile Fine Tip w/ Ruler & 9 Blank Labels	100

# UTILITY MARKERS



## 1453 WATERPROOF PERMANENT MARKING MINI MARKER

Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1453-100	Non-Sterile Regular Tip	100
1453SRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



## 1462 DUAL TIP BARREL XL PREP RESISTANT INK SKIN MARKER


Regular Tip + XL Prep Resistant Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1462XLSR-100	Sterile w/ Ruler 1 Regular & 1 Fine Tip	100
1462XLSRL9-100	Sterile 1 Regular & 1 Fine Tip w/ Ruler & 9 Blank Labels	100



## PRE-PRINTED MEDICATION LABELS

Custom labels available for no additional cost.



# PRE-PRINTED MEDICATION LABELS

Ancef	Lidocaine _____ %
Ancef	Lidocaine _____ %
Adenosine _____ mcg/ml	Nipride _____ mcg/ml
Adenosine _____ mcg/ml	Nitro _____ mcg/ml
Contrast Name _____ mgI/ml	Normal Saline 0.9%
Contrast Name _____ mgI/ml	Normal Saline 0.9%
Heparin _____ units/ml	Verapamil _____ mcg/ml
Heparin _____ units/ml	Verapamil _____ mcg/ml
Heparinized Saline _____ units/ml	
Heparinized Saline _____ units/ml	

## CATH LAB COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
VCATH-20	Sterile 20 Color Kit	100
VCATH-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

AIR	Mitomycin _____ mg/ mL
Antibiotic	Moxifloxacin 0.5% (Vigamox)
BSS	Perfluoron
BSS	SF6 _____ %
C3F8 _____ %	Sterile Water
Decadron _____ mg/mL	Triamcinolone PF _____ mg/mL
Dexamethasone _____ mg/mL	Trypan Blue
Lidocaine _____ % MPF	Vitrase _____ units/mL
Miochol	
Miostat	

## OPHTHALMIC COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
VEYE-20	Sterile 20 Color Kit	100
VEYE-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

# PRE-PRINTED MEDICATION LABELS

<b>Contrast</b> Name _____ mgI/ml	<b>Bupivacaine</b> _____ %
<b>Contrast</b> Name _____ mgI/ml	<b>Bupivacaine</b> _____ %
<b>Heparin</b> _____ units/ml	<b>Lidocaine 1% w/ Sodium Bicarb</b> Amt _____ Exp _____
<b>Heparin</b> _____ units/ml	<b>Lidocaine 1% w/ Sodium Bicarb</b> Amt _____ Exp _____
<b>Heparinized Saline</b> _____ units/ml	<b>Normal Saline</b> 0.9%
<b>Heparinized Saline</b> _____ units/ml	<b>Normal Saline</b> 0.9%
<b>Lidocaine</b> _____ %	<b>1% Lidocaine HCL w/ Epinephrine 1:100,000</b> Amt _____ Exp _____
<b>Lidocaine</b> _____ %	<b>1% Lidocaine HCL w/ Epinephrine 1:100,000</b> Amt _____ Exp _____
<b>Nipride</b> _____ mcg/ml	
<b>Nitro</b> _____ mcg/ml	

## INTERVENTIONAL RADIOLOGY COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
<b>VIR-20</b>	Sterile 20 Color Kit	100
<b>VIR-22R-20</b>	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

<b>Bupivacaine</b> _____ %	<b>Depomedrol</b> _____ mg/ml
<b>Bupivacaine</b> _____ %	<b>Depomedrol</b> _____ mg/ml
<b>Lidocaine</b> _____ %	<b>Betamethasone</b> _____ mg/ml
<b>Lidocaine</b> _____ %	<b>Betamethasone</b> _____ mg/ml
<b>Contrast</b> Name _____ mgI/ml	<b>Dexamethasone</b> _____ mg/ml
<b>Contrast</b> Name _____ mgI/ml	<b>Dexamethasone</b> _____ mg/ml
<b>Kenalog</b> _____ mg/ml	<b>Verapamil</b> _____ mcg/ml
<b>Normal Saline</b> 0.9%	<b>Verapamil</b> _____ mcg/ml
<b>Normal Saline</b> 0.9%	
<b>Triamcinolone</b> _____ mg/mL	

## PAIN MANAGEMENT COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
<b>VPAIN-20</b>	Sterile 20 Color Kit	100
<b>VPAIN-22R-20</b>	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

# PRE-PRINTED MEDICATION LABELS

<b>Acetic Acid</b>	<b>Methergine</b>
<b>Bupivacaine _____% with Epinephrine</b>	<b>Methylene Blue 50mg/5ml</b>
<b>Bupivacaine _____%</b>	<b>Mineral Oil</b>
<b>Estradiol 0.01% Cream</b>	<b>Monseal's Solution</b>
<b>Hydrogen Peroxide 3%</b>	<b>Normal Saline 0.9%</b>
<b>Iodine Scrub</b>	<b>Silver Sulfadiazine 1% (Silvadene)</b>
<b>Iodine Prep</b>	<b>Vasopressin 20u/ml</b>
<b>KY Jelly</b>	<b>WATER</b>
<b>Lidocaine 1% Plain</b>	
<b>Lugol's Solution</b>	

## OBGYN COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
<b>VOBGYN-20</b>	Sterile 20 Color Kit	100
<b>VOBGYN-22R-20</b>	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

<b>Antibiotic</b> Amt _____ Exp _____	<b>Lidocaine _____%</b> Amt _____ Exp _____
<b>Bacitracin</b> Amt _____ Exp _____	<b>Marcaine _____%</b> Amt _____ Exp _____
<b>Cocaine _____%</b> Amt _____ Exp _____	<b>Marcaine _____% with Epinephrine</b> Amt _____ Exp _____
<b>DYE</b> Amt _____ Exp _____	<b>Lidocaine _____% with Epinephrine</b> Amt _____ Exp _____
<b>Epinephrine</b> Amt _____ Exp _____	<b>Peroxide</b> Amt _____ Exp _____
<b>Heparin _____units</b> Amt _____ Exp _____	<b>Normal Saline 0.9%</b> Amt _____ Exp _____
<b>Heparinized Saline units/ml</b> Amt _____ Exp _____	<b>Water</b> Amt _____ Exp _____
<b>Miochol</b> Amt _____ Exp _____	<b>Bupivacaine _____%</b> Amt _____ Exp _____
(BLANK)	(BLANK)
(BLANK)	(BLANK)

## GENERAL O.R. LABELS (COLOR)

PART #	DESCRIPTION	CASE QTY
<b>V VSC04CL</b>	Sterile 20 Color Kit	100
<b>V20063DC-DTR</b>	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100
<b>V20063DC</b>	Sterile 20 Color Kit w/ 1408 Reg Tip Permanent Ink Marker	100
<b>V20063DC-F</b>	Sterile 20 Color Kit w/ 1411 Fine Tip Permanent Ink Marker	100
<b>V20063DC-S</b>	Sterile 20 Color Kit w/ 1416 Ultra Fine Tip Permanent Ink Marker	100

# PRE-PRINTED MEDICATION LABELS

<b>Antibiotic</b> Amt _____ Exp _____	<b>Lidocaine _____%</b> Amt _____ Exp _____
<b>Bacitracin</b> Amt _____ Exp _____	<b>Marcaine _____%</b> Amt _____ Exp _____
<b>Cocaine _____%</b> Amt _____ Exp _____	<b>Marcaine _____% with Epinephrine</b> Amt _____ Exp _____
<b>DYE</b> Amt _____ Exp _____	<b>Lidocaine _____% with Epinephrine</b> Amt _____ Exp _____
<b>Epinephrine</b> Amt _____ Exp _____	<b>Peroxide</b> Amt _____ Exp _____
<b>Heparin _____units</b> Amt _____ Exp _____	<b>Normal Saline 0.9%</b> Amt _____ Exp _____
<b>Heparinized Saline _____units/ml</b> Amt _____ Exp _____	<b>Water</b> Amt _____ Exp _____
<b>Miochol</b> Amt _____ Exp _____	<b>Bupivacaine _____%</b> Amt _____ Exp _____
(BLANK)	(BLANK)
(BLANK)	(BLANK)

## GENERAL O.R. LABELS (BLACK & WHITE)

PART #	DESCRIPTION	CASE QTY
<b>V20063D</b>	Sterile 20 B&W Kit w/ 1408 Permanent Ink Marker	100



## BLANK MEDICATION LABELS

PART #	DESCRIPTION	CASE QTY
<b>V20063DB-DTR</b>	Sterile 10 Blank Labels w/ 1422 Dual Tip Marker & Ruler	100

# PRE-PRINTED MEDICATION LABELS

## EXTRA VESSELS

FOR USE IN ORGAN TRANSPLANTATION ONLY



DONOR ID

DONOR ABO

RECOVERY DATE

### INFECTIOUS DISEASE TEST RESULTS

Extra vessels are from a donor with:

1. positive results (NAT included) for HIV, HBV (HBsAg) or HCV\*

YES  NO  PENDING

2. positive results for total anti-HBc

YES  NO  PENDING

3. HIV, HBV, HCV Criteria

YES  NO

\*OPTN Policy prohibits storage of these extra vessels. They must be quarantined and destroyed if not used in intended recipient.

Check source documents for results.

11/2020

OPTN ORGAN PROCUREMENT AND TRANSPLANTATION NETWORK

## VESSEL LABELS

### PART #

### DESCRIPTION

### CASE QTY

**V21VESSEL-2X4**

2" x 4" Sterile Vessel Labels

100

**V21VESSEL-16-2X4**

2" x 4" Sterile Vessel Labels w/ 1416 Ultra Fine Tip Permanent Ink Marker

100

**V21VESSEL-2X5**

2" x 5" Sterile Vessel Labels


100

**V21VESSEL-16-2X5**

2" x 5" Sterile Vessel Labels w/ 1416 Ultra Fine Tip Permanent Ink Marker

100

## 1" diameter Circle



## KIDNEY TAGS

### PART #

### DESCRIPTION

### CASE QTY

**9PHKND96W**

Sterile 1" Round Polypropylene Tags, Left (Yellow) & Right (Blue) Kidneys

100

**9PHKND105W**

Sterile 1" Round Polypropylene Tag, Left Kidney (Yellow)

100

**9PHKND106W**

Sterile 1" Round Polypropylene Tag, Right Kidney (Blue)

100

# CUSTOM MEDICATION LABELS

Black & White or Color Coded

Lidocaine ___%	<b>Contrast</b>	Heparin ___ units 0.9% Saline ___ ml
Lidocaine ___%	Adenosine ___ mcg/ml	0.9% Saline ___ cc
Lidocaine ___%	Adenosine ___ mcg/ml	0.9% Saline ___ cc
Lidocaine ___%		0.9% Saline ___ cc
Lidocaine ___%		0.9% Saline ___ cc
Lidocaine ___%		Ancef ___ grams + 0.9% Saline
Marcaine ___%		Ancef ___ grams + 0.9% Saline
Marcaine ___%		Marcaine ___ grams + 0.9% Saline
Marcaine ___%	Heparin ___ units 0.9% Saline ___ ml	Marcaine ___ grams + 0.9% Saline
<b>Contrast</b>	Heparin ___ units 0.9% Saline ___ ml	

Epinephrine 1:100,000 By Use Date/Time _____	Marcaine 0.25% By Use Date/Time _____
Epinephrine 1:100,000 By Use Date/Time _____	Marcaine 0.25% By Use Date/Time _____
Epinephrine 1:200,000 By Use Date/Time _____	Marcaine 0.5% By Use Date/Time _____
Epinephrine 1:200,000 By Use Date/Time _____	Marcaine 0.5% By Use Date/Time _____
Lidocaine 0.5% By Use Date/Time _____	Mepivacaine 1.5% By Use Date/Time _____
Lidocaine 0.5% By Use Date/Time _____	Mepivacaine 1.5% By Use Date/Time _____
Lidocaine 1% By Use Date/Time _____	Normal Saline 0.9% By Use Date/Time _____
Lidocaine 1% By Use Date/Time _____	Ropivacaine Strength: _____ By Use Date/Time _____
Lidocaine 2% By Use Date/Time _____	Ropivacaine Strength: _____ By Use Date/Time _____
Lidocaine 2% By Use Date/Time _____	Ropivacaine Strength: _____ By Use Date/Time _____

Bupivacaine 0.5%
Bupivacaine 0.25%
Depomedrol 80mg/ml
Dexamethasone 4mg/ml
Omnipaque 300
Omnipaque 240
Lidocaine 1%
Lidocaine 2%
Normal saline

## Waterproof & fluid resistant

- Ink will not come off or become illegible due to contact with fluids

## Smudge Proof

- Eliminates concerns of errors due to smearing

## Residue Free

- Eliminates build-up on stainless basins and cups

## Tear Free

- Simplifies cleanup for sterile processing, saving time

**VLSCOT**<sup>®</sup>  
MEDICAL, LLC.

## **SHARPS & STAFF SAFETY**

Reusable - Autoclavable



# SOFFZONE® NEUTRAL ZONE\*



**3" x 9" SOFFZONE® NEUTRAL ZONE MAT**

PART #	DESCRIPTION	CASE QTY
<b>FNZ309RED</b>	3"x9" Red Tray w/ Raised Contours for Sharps	6



**6" x 9" SOFFZONE® NEUTRAL ZONE MAT**

PART #	DESCRIPTION	CASE QTY
<b>FNZ609RED</b>	6"x9" Red Mat w/ Raised Contours for Sharps	6



**SOFFZONE® SYRINGE TRAY**

PART #	DESCRIPTION	CASE QTY
<b>FSY309RED</b>	3"x9" Red Tray w/ Tabs for Holding Syringes	6

\*Soffzone® products are covered by US Patents 10111710 B2, 10905515 B2 and other patents pending.



## **ADDITIONAL PRODUCTS**



# URINALS



**MILLIE FEMALE URINAL**

PART #	DESCRIPTION	CASE QTY
201000-12	Female Urinal	12
201036	Female Urinal	36



**UNIVERSAL URINAL w/ LID**


PART #	DESCRIPTION	CASE QTY
201000-L-12	Universal Urinal w/ Leak Resistant Lid - Individually Bagged w/ Instructions	12
201036-L	Universal Urinal w/ Leak Resistant Lid - Individually Bagged w/ Instructions	36

# RULERS

## Rulers & labels included in sterile packaging:

Viscot Medical Sterile markers come individually packaged often with a fluid resistant ruler and sometimes with back up labels. Rulers are flexible, waterproof, and include both metric and U.S. measurements units (15 cm / 6 in).

Blank back up labels included with markers come in quantities of 4, 6 or 9. They are paper based labels and are designed as a back to our standard pre-printed medication labels.



### 1410 METRIC/INCH (15CM/6IN) RULER

PART #	DESCRIPTION	CASE QTY
1410-100	Non-Sterile	100
1410-1000	Non-Sterile	1,000
1410S-300	Sterile	300

### R24 24 INCH FOLDED RULER

PART #	DESCRIPTION	CASE QTY
R24-1000	Non-Sterile Paper Ruler	1,000

# DRAPES

Viscot Medical carries a selection of non-sterile procedural drapes that come Poly Lined or Poly backed.

- Available Material: Tissue Poly Tissue, Plastic, and Nonwoven
- Available Sizes: 17x19, 18x26, 24x26 & 36x24
- Options of a 3" and 4" round fenestration available with or without adhesive patch/ tape



## STERILE MINOR SURGICAL DRAPES









PART #	DESCRIPTION	CASE QTY
1503	18" x 26" Poly-Lined White Drape	300
1501	18" x 26" Poly-Lined White Drape, Fenestrated - 2 3/4"	300

## CSR WRAPS



- Available Material:  
2 layer Cellulose base 260 grade material
- Available Sizes:  
18", 20", 24", 30", 36", 40" and 48"

# HOW TO READ YOUR PACKAGING

WHAT IS THIS SYMBOL?	WHAT IT MEANS:	WHERE WILL I SEE IT?
	General warning sign	To signify a general warning
	CE marking	A device is in conformity with the applicable requirements set out in this Regulation and other applicable Union harmonisation legislation providing for its affixing
	Do not use if package is damaged and consult instructions for use	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information
	Authorized Representative in the European Community/ European Union	Indicates the authorized representative in the European Community / European Union
	Use by date	Indicates the date after which the medical device is not to be used
	Instructions for Use	Indicates that there are instructions available for this product
	Latex Free	Indicates the item is Latex Free
	Manufacturer	Indicates the medical device manufacturer

# HOW TO READ YOUR PACKAGING

	WHAT IS THIS SYMBOL?	WHAT IT MEANS:	WHERE WILL I SEE IT?
	Medical device	Indicates the item is a medical device	ISO/DIS 15223- 1:2021 Reference no. 5.7.7
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process	ISO 15223-1: 2021 Reference no. 5.2.7. (ISO 7000-2609)
	Catalogue number Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified	ISO 15223-1: 2021 Reference no. 5.1.6. (ISO 7000-2493)
	Batch code/ Lot number	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1 Reference no. 5.1.5
	Do not re-use	Indicates a medical device that is intended for one single use only NOTE: Synonyms for "Do not reuse" are "single use" and "use only once"	ISO 15223-1:2021 Reference no. 5.4.2. (ISO 7000- 1051)
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1:2021 Reference no. 5.2.4. (ISO 7000-2502)
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.	ISO 15223-1:2021 Reference no. 5.2.3. (ISO 7000-2501)
	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed	ISO 15223-1: 2021 Reference no. 5.3.6. (ISO 7000-0533)



[www.viscot.com](http://www.viscot.com)