

Viper Plus Reclining Wheelchair

HCPCS: K0004, E1226-1, K0004, E1226-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Patented built-in seat rail extension and extendable upholstery easily adjusts seat depth from 16" or 18"
- Attendant wheel locks in a convenient and easily accessible location
- Spreader bar acts as a back stabilizer and a convenient push handle for the attendant
- Offset rear wheel position adds additional stability and equally distributes the patient's weight while in the various reclining positions
- Hydraulic reclining mechanism allows for adjustments up to 180°
- Durable, flame-retardant nylon upholstery resists mildew and bacteria
- Carbon steel frame with double-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels are lightweight and maintenance free also features 8" front casters
- Nylon, padded armrests provide added user comfort, and feature the flip-back style for easy transfers
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Standard dual axle provides easy adjustment of seat-to-floor height to hemi-level
- Front caster forks are adjustable in three positions
- Comes standard with swing-away elevating leg rests
- Composite foot plates are durable and lightweight
- Polyurethane rear tires and casters offer superior performance, smoother ride and are lightweight

PRODUCTS

Description

PLA416RBDDA	ViperGT,16"Recl,FlpBackRemDsk,1/ea
PLA416RBDFA	ViperGT,16"Recl,FlpBackRemFull,1/ea
PLA418RBDDA	ViperGT,18"Recl,FlpBackRemDesk,1/ea
PLA418RBDFA	ViperGT,18"Recl,FlpBackRemFull,1/ea
PLA420RBDDA	ViperGT,20"Recl,FlpBackRemDesk,1/ea
PLA420RBDFA	ViperGT,20"Recl,FlpBackRemFull,1/ea

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Product Weight	47 lb 49 lb 51 lb
Product Dimensions	42" (L) w/ riggings x 24" 26.5" 28" (W) x 50" (H)
Seat Dimensions	16" 18" 20" (W) x 16" - 18" (D)
Weight Capacity	300 lb
Seat To Floor Height	17.5" 19.5"
Carton Shipping Weight	53 lb 55 lb 57 lb
Folded Dimensions	42" (L) w/ riggings x 12.5"(W) x 50" (H)
Back Height	23" 33" w/ headrest
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov]

ACCESSORIES

Description

STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX
WASL	Limb Support,Left,1/cs, CS
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA