



VIVE 8" ALTERNATING PRESSURE MATTRESS

LVA1066

HCPCS: E0277

Designed to be used in the prevention, treatment, and management of pressure ulcers, spots, and bedsores, the alternating pressure mattress system includes a 16-cell mattress and variable pressure pump. Perfect for bedridden patients suffering from bedsores, pressure spots, and ulcers.

SPECIFICATIONS

Dimensions:

- Length: 80"
- Width: 34"
- Height: 8"

What's included:

- Vive Alternating pressure mattress
- Vive digital compressor pump
- Connecting air hose
- Replacement air chamber
- Patching Kit
- One year guarantee on pump

Care Instructions:

- For best results, wipe down with a damp cloth coated with a detergent product solution or disinfectant.
- Do not use corrosive cleaning products such as industrial degreasers or acetic solvents.

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ALLEVIATES BED SORES AND ULCERS

Alleviates pressure points while relieving ulcers.



SIMPLE TOUCH CONTROLS

Customize air flow, firmness and positioning.



ULTRA QUIET PUMP

Provides restful sleep and hooks onto bed.

PREVENT BEDSORES AND ULCERS

Alleviates painful pressure on patients who have limited mobility or are bedridden.

19 ALTERNATING AIR CHAMBERS

Continually redistributing the weight, the air chambers inflate and deflate alternatively.

VARIABLE PRESSURE SETTINGS

Adjust the amount of air pressure for personal comfort and support or select the static or seated setting as necessary.

ULTRA QUIET COMPRESSOR

Easy to set up, the digital compressor pump is exceptionally quiet to provide a restful atmosphere.

REMOVABLE WATERPROOF COVER

Constructed with a durable, latex-free vinyl, the waterproof alternating mattress is easy to clean.

8" TWIN XL MATTRESS

Supporting up to 350 pounds, the mattress is engineered for even support and breathable comfort.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - BR3