

vive health

portal.vivehealth.com



VIVE

TRANSPORT STAND ASSIST MOB1043WHT

Ideal for aiding those who can support a portion of their own weight, the Vive transfer stand assist allows individuals to stand and be transported to a bed, wheelchair, toilet, or chair without the use of a transport sling. Minimizing the risk of hazardous falls, the transfer stand assist is generously padded for patient and caregiver comfort. The durable steel frame also has many safety features including dual locking casters, locking swivel seat pads, a wide footrest, and support strap hooks for caregiver peace of mind.

WHAT'S INCLUDED

- 1pc Transport Stand Assist

**IDEAL FOR AIDING
THOSE WHO CAN SUPPORT A
PORTION OF THEIR OWN WEIGHT**



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com

Specifications

Dimensions	44" h x 30" w x 25" l
Footrest base	16" x 12"
Inner Diameter for Patient Space	19.75" w x 18" l
Seat Size	19"
Weight	60 lb
Material(s)	Steel, Polymer, and Foam.
Weight limit	400LBS



AID IN TRANSFERS, STANDING, AND LIFTING

Designed for easy use by individuals and caregivers, the Vive transfer stand assist provides mobility assistance for those who can support a portion of their own body weight. The transfer stand assist minimizes the risk of falls and caregiver injury at home or in skilled nursing facilities when transferring from a bed to a wheelchair, toilet, sofa, and more.

PADDED FOR GREATER COMFORT

Each handrail is enclosed in soft foam to provide a secure, yet comfortable grip that minimizes fatigue. The shin rests and split seat pads are also cushioned for patient comfort.

DURABLE STEEL FRAME

The Vive transfer stand assist is constructed with a strong steel frame that is perfectly sized to fit through most doorways to fit through most doorways, and safely supports up to 400 pounds. The transfer stand assist rests on four swivel casters that easily glide across multiple flooring surfaces including laminate and carpeting.

INTEGRATED SAFETY FEATURES

Securely locking in place, the transfer assist's dual locking brakes eliminate unwanted movement while standing or transporting a patient or loved one. The frame also includes hooks to easily attach side support straps for added security and stability.

VIVE GUARANTEE

A two-year guarantee so you can purchase now with confidence.

