

VIVE

CORE WHEELCHAIR

MOB1100

The Core Wheelchair redefines mobility with its perfect balance of strength, comfort, and convenience. Weighing just 39 pounds, its lightweight yet durable steel construction is designed to empower your active lifestyle. With an ergonomic 18" x 16" seat and an ideal seat height of 19.5", this wheelchair offers exceptional comfort and support for users of all sizes.

Engineered for versatility, the Core Wheelchair features 24-inch rear wheels for smooth performance on a variety of terrains. Its compact, travel-friendly design folds to an ultra-slim 11 inches, making storage and transportation effortless. Whether you're navigating daily errands or heading out for an adventure, the Core Wheelchair is the perfect companion.

SPECIFICATIONS

What's Included:

➤ Core wheelchair

Color(s):

➤ Black or Blue

Product Dimensions:

➤ Seat: 18" W x 16" D

➤ Seat Height: 19.5"

➤ Overall dims: 43L x 36H x 25"W

➤ Folded width: 11"

➤ Front wheels: 8"

➤ Rear wheels: 24"

Product Weight:

➤ 39lbs

Product Material(s): Steel



Black



Blue

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



DURABLE FOLDING FRAME

Built with high-quality steel, this wheelchair combines durability and lightweight portability. Weighing 39 lbs, it folds down to just 11" in width, making storage and transport hassle-free.



OPTIMIZED COMFORT DESIGN

With an 18" W x 16" D seat and a 19.5" seat height, this wheelchair ensures a supportive and comfortable experience for standard-sized users.



RELIABLE CONSTRUCTION AND WHEELS

Equipped with 8" front wheels and sturdy 24" rear wheels, this chair glides smoothly over various surfaces, offering stability and control.

SAFE AND CONVENIENT FEATURES

This latex-free wheelchair is designed with safety in mind, providing users with a secure and comfortable ride for everyday use.

CUSTOMIZABLE COMFORT

Adjustable legrests and swingback armrests for a personalized, flexible seating experience.



(239) 220-5772



sales@vivehealth.com

vive
health

SR25-1R1