



## VIVE

### NEPTUNE MOBILITY SCOOTER

MOB1126MBL

The Vive Neptune Mobility Scooter is engineered for smooth travel across a variety of indoor and outdoor surfaces. With 9" flat-free tires, 2.5" of ground clearance, and anti-tip rear wheels, the Neptune glides over door thresholds, sidewalks, and uneven terrain with ease. A powerful quick-connect battery system delivers up to 16 miles of range on a single charge, while the 3.75 mph top speed ensures steady, reliable performance. Designed with a wide, padded seat, ergonomic backrest, and adjustable armrests, the Neptune offers customizable comfort and supports users up to 300 pounds. The user-friendly delta tiller includes simple speed controls, a bright LED headlight, and an easy-to-read battery indicator. For added convenience, the scooter disassembles into five lightweight pieces for transport, including the front basket to keep essentials within reach.

#### WHAT'S INCLUDED

- 1x Neptune Mobility Scooter.
- 1x Storage Basket.

**SAFELY AND SMOOTHLY  
NAVIGATING OVER ANY SURFACE**



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com

#### GENERAL

Feature	Value
Maximum weight capacity	300 lbs (136 kg)
Drive type	Rear wheel drive
Maximum built-in speed limit	3.75 mph (6 km/h)
Turning radius	62.99 inches (160.00 cm)
Ground clearance	2.5 inches (6.3 cm)
Maximum Travel Distance	16 Miles
Brake type	Electro-mechanical

#### SEAT

Feature	Value
Seat Type	Stadium style (seat-back foldable)
Seat Size (depth x width)	15 x 17.5 inches (38.1 x 44.4 cm)
Seat Color	Black
# of seat height positions	4
Height from deck to lowest seat position	15 inches (38.1 cm)
Height from deck to highest seat position	18 inches (45.7 cm)

#### MOTOR

Feature	Value
Product weight	116 lbs (52.6 kg)
Motor type	DC
Total number of motors	1
Voltage	24V
Wattage	180W
Control type	Dual throttle control levers
Ignition key	Console

#### BATTERY

Feature	Value
Battery pack type	Lead-acid
Total # of batteries	2 (in 1 sealed battery pack housing)
Removeable?	Yes
Battery weight	29.5 lbs (13.3 kg)
Battery voltage	2 x 12V
Battery amp-hours	2 x 20Ah
Total charging time	Up to 8 hours
Total battery lifespan	Approx. 300 charging cycles
Type of charging	Offboard
Level indicator	Yes (on console)
Charging indicator	Red = in progress, Green = charge complete

#### WHEELS

Feature	Value
Front wheel diameter	8 inches (20.3 cm)
Rear tire/wheel diameter	9 inches (22.8 cm)
Anti-tip wheel diameter	3 inches (7.6 cm)
Type of tire/wheel	Solid
Total # of tires/wheels	4 (+ 2 anti-tip wheels)
Maximum incline limit ***	9 °

## EFFORTLESS ALL-TERRAIN MOBILITY

Navigate smoothly over carpet, tile, sidewalks, grass, and uneven paths with 9" flat-free tires and 2.5" ground clearance. Rear anti-tip wheels provide added control on inclines, making indoor and outdoor travel safer and more stable.

## EXTENDED BATTERY LIFE

Ride up to 16 miles per charge with a quick-connect battery that's easy to remove for fast charging. A top speed of 3.75 mph ensures steady, reliable performance for daily errands and outings.

## CUSTOMIZABLE FOR GREATER COMFORT

Sit back in a wide, cushioned seat with a supportive backrest and height-adjustable design. The armrests flip up for easy transfers and can be positioned for a personalized, ergonomic fit.

## INTUITIVE STEERING

The delta tiller is intuitive and easy to operate, even with limited hand mobility. Monitor speed and battery level at a glance and ride safely with the built-in LED headlight for low-light visibility.

## DISASSEMBLES FOR EASY TRAVEL

Quickly break down the scooter into five manageable pieces for car transport or compact storage. A detachable front basket offers added convenience for shopping, travel, or daily use.

