

SKU: RHB3003

coretech®
by vivehealth

Rx

LYMPHEDEMA PUMP



The Coretech RHB3003 Lymphedema Pump is a pneumatic compression therapy device specifically designed to address lymphedema and related vascular conditions. This 510K-cleared medical device utilizes sequential compression technology to promote lymphatic drainage and improve circulation in patients suffering from lymphatic impairment. By delivering precise, graduated pressure patterns, the system helps redirect excess lymphatic fluid from compromised areas to healthy regions where it can be naturally processed by the body, providing patients with effective symptom relief and improved quality of life.

20-80 MMHG 5 MMHG INCREMENTS



8-CHAMBERS SLEEVE DESIGN

8-Chamber sleeve configuration with 4 air inlets and shared chambers create optimal segmental air pressure therapy for enhanced lymphatic drainage and superior patient outcomes.



BILATERAL AND UNILATERAL TREATMENT

Included air adapter enables flexible treatment configurations - treat one leg, both legs simultaneously, or combine leg and arm therapy to address multiple affected areas in a single session.



PATIENT SAFETY LOCK FUNCTION

Built-in locking mechanism ensures patient safety and preventing therapy interruption, providing uninterrupted treatment sessions for optimal therapeutic results



DUAL MODE THERAPY OPTIONS

Choose between Peristaltic Mode (Mode A) and Sequential Gradient Mode (Mode B). Each mode produces distinct compression patterns optimized for different patient needs and treatment protocols. Select from 30, 45, and 60-minute treatment times.

SKU: RHB3003

WHAT'S INCLUDED:

- 1 RHB3003 Pump
- 2x Leg Sleeves (8 Chambers)
- 2x Sleeve Extenders
- 2x Air Tubes
- 1X Air Adapter
- 1X Allen Key



SIZING GUIDE – LEG SLEEVES

	CIRCUMFERENCES					
Size	Thigh Below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch
Small	31"	30"	28"	24.5"	21.5"	19.5"
Medium	32"	30"	28"	24.5"	21.5"	19.5"
Medium Wide	36.5"	34"	30"	28"	25.5"	22.5"
Large	33"	32"	30"	26"	21.5"	19.5"

PRODUCT MATERIAL(S):

Cuff: Nylon TPU (70/30)

Pump: Plastic/Metal (50/50)

LATEX FREE: Yes

CARE: Spot clean with a damp cloth and dry immediately

