
CENTURY **GOLD**

The **Century Gold** mattress system provides alternating pressure therapy with true low air loss to aid patients who are at risk of developing pressure ulcers or have existing stage 1, 2, 3, or 4 ulcers. Caregiver friendly, the Century Gold is an affordable, therapeutic surface solution for any clinical care setting.







On demand raised side rail feature for added patient safety.

MASTER CONTROL UNIT	
Dimension:	(L)15.5"x (W)7.5"x (H)9.5"
Weight:	14.5 lbs
Rating Voltage:	110-240 volts
Max Current:	5A max.
Rating Frequency:	60Hz
Cycle Time:	10, 15, 20, 30 mins
	Static function
	2 in 1 alternation
Min/Max Pressure:	15~50 mmHg +/- 6 mmHg
Air flow rate:	> 30L/min
Protection Type:	Class I, Type BF
Low Pressure Alarm:	Visual / Audio
Power Failure Alarm:	Visual / Audio

MATTRESS REPLACEMENT	
Model No:	Century Gold
Dimensions:	(L)80"x (W)36", 42", 48" x (H)10"
Weight:	25.3 lbs
Weight Capacity:	500 lbs (W)36", 700 lbs (W)42", 750lbs (W)48"
Cell Material:	Clear TPU
Cell Number:	21 cells
Cover Material:	Blue Nylon with PU coating finish
Base Material:	Black Polyester with PVC coating
	CPR Quick Release System
	Cable Management System on bottom cover
HCPSC Code	E0277

PRIUS HEALTHCARE USA
Creating innovation and difference

T (813) 854-5464
E info@priushc.com
W priushcusa.com

160 SCARLET BLVD
OLDSMAR, FL 34677

CENTURY **GOLD**

Provides you, the caregiver, with peace of mind, and your patient with world class care and comfort



PRIUS HEALTHCARE USA
Creating innovation and difference

Perfect for patients and caregivers




On demand raised side rail for added patient safety



On demand true low air loss for moisture management



Slimline pump design is unobstructive



PRESSURE MAPPING

Century Gold provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.



CENTURY GOLD

- Interchangeable air cell design provides easy cleaning and minimizes the proliferation of bacteria
- The Century Gold employs 8" air bladders to provide therapeutic 2-1 alternating pressure and On Demand True Low Air Loss
- The Century Gold replacement mattress can accommodate patient weights up to 750 lbs
    - During a power outage the system stays closed allowing the air cells to remain inflated
        - Air-filled safety side rails



Century Gold is your partner in prevention and treatment of pressure ulcers



Cable management system keeps electrical cables off floor for caregiver and patient safety



CPR quick release system deflates the system within 20 seconds



Spring loaded hooks snap to fit your bed's footboard

Prius