

DYNAFLOW IMMERSION

DynaFlow Mattress System is a low to medium risk system for patients who are at risk of developing Pressure Ulcers according to your sound clinical judgment. DynaFlow can also be used for patients who have an existing stage 1, 2, 3 and 4 pressure ulcers, in conjunction with your hospitals policy on pressure area management.





MASTER CONTROL UNIT	
Dimension:	(L)12.3"x (W)6.7"x (H)4.9"
Weight:	14.9 lbs
Rating Voltage:	110-120 volts
Max Current:	5A
Rating Frequency:	60Hz
Cycle Time:	3~95 mins
	Static / Pulsation / Immersion function
	3 in 1 alternation
Min/Max Pressure:	15~52 mmHg +/- 4 mmHg
	10 Digital Adjustable Comfort Control
Protection Type:	Class I, Type BF
Bottom Out Alarm:	Visual / Audio
Power Failure Alarm:	Visual / Audio
Features:	User-Friendly LCD Control Panel
	Fowler Function
	Auto-Set Function

MATTRESS REPLACEMENT	
Model No:	DynaFlow
Dimensions:	(L)80"x (W)36", 42", 48" x (H)10"
Weight:	39.6 lbs
Weight Capacity:	500, 750, 1000 lbs
Cell Material:	Black Nylon with PU Coating
Cell Number:	18 cells
Cover Material:	Blue Nylon with PU coating finish
Base Material:	Black Polyester with PVC coating
	CPR Quick Release System

PRIUS HEALTHCARE USA
Creating innovation and difference

T (813) 854-5464
E info@priushc.com
W priushcusa.com

160 SCARLET BLVD
OLDSMAR, FL 34677



DYNAFLOW
IMMERSION

Premium Low Air Loss and Alternating
Pressure Mattress Replacement System



PRIUS HEALTHCARE USA
Creating innovation and difference

priushcusa.com


Perfect for patients and caregivers




Immersion mode provides the optimal pressure-reducing therapy



The modular construction of the cells makes for simplistic repair and cleaning



Incorporated with automatic sensing technology for fowler detection and pressure adjustment



PRESSURE MAPPING

DynaFlow provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.

DYNAFLOW
IMMERSION




- Low Air Loss, Pulsation, 3-1 Alternation and Immersion Therapies in one system
- Nylon PU cell design provides easy cleaning and minimizing the proliferation of bacteria
- Mattress replacement can accommodate patient weights up to 1,000 lbs
- Equip with 2" convoluted foam base provides additional support during power failure situation

Automatically sets optimal pressure based on patient's height and weight.

60 seconds of Pulsation simulates a massage to assist in maximizing a patients comfort.

1-5 minutes of Immersion therapy



DynaFlow can detect changes in the angle of the bed frame with "Fowler Mode" and adjust the pressure based on these results

Incorporates 4 modes of Low Air Loss therapy: Static, Pulsation, 3-1 Alternation and Immersion.

Prius

DynaFlow is your best choice for protection and treatment of pressure ulcers



Cable management keeps electrical cables off the floor



CPR quick release system deflates the system within 20 seconds



Hanging bracket fits on all type of beds