PRIUS ULTRA CONVERTIBLE



THE PRIUS ULTRA CONVERTIBLE IS A SELF-ADJUSTING PRESSURE RELIEVING MATTRESS SYSTEM. INTERNAL ZONED MULTI LAYERED FOAM AIR CELLS USE INTEGRATED CHAMBER TECHNOLOGY CONSTANTLY MIGRATING AIR THROUGHOUT THE SYSTEM WHILE SELF-ADJUSTING TO EVERY PATIENT'S WEIGHT, SHAPE AND SIZE.

FEATURES

- SERENE® FOAM MATTRESS TOPPER IS THE NEWEST MOST ADVANCED FOAM AVAILABLE PROVIDING UNPRECEDENTED COMFORT AND MAXIMUM IMMERSION TO THE MOST HIGH-RISK PATIENTS.
- MULTI-LAYER FOAM HEEL SECTION OFFLOADS LOWER LEG PRESSURE TO FLOAT THE SENSITIVE HEEL AREA.
- FULL LENGTH LOWER RAIL BOLSTERS SUPPORTS KEEP PATIENT CENTERED ON MATTRESS SURFACE WHILE ASSISTING IN BED INGRESS AND EGRESS.
- OPTIONAL RAISED BOLSTERS AVAILABLE.
- TOP COVER FABRIC IS A MICRO-CLIMATE, VAPOR PERMEABLE, FLUID PROOF, LOW FRICTION 4-WAY STRETCH IMMERSION FABRIC THAT PROMOTES LOW FRICTION AND SHEARING.
- OPTIONAL POWERED ALTERNATING PRESSURE CONTROL UNIT FOR ENHANCED TREATMENT AND WOUND HEALING.

SPECIFICATIONS

MATTRESS ZONE SYSTEM: head, shoulder, torso, foot
WIDTH: 36″, 39", 42″, or 48"
LENGTH: 80″ or 84″
HEIGHT: 7″
WEIGHT CAPACITY: 600 lbs to 1,000 lbs
FLAMMABILITY: Federal Standard 16 CFR 1633
WARRANTY: Pump - 1 Year, Mattress - 5 Years

ITEM NUMBER

SYSTEM - PHR-CHM36
MATTRESS - FM-PHRCHM3680
PUMP - FC-PHR0052
HOSE - PM-PHR0335

